



Date: September 15, 2022

Author: Software Engineering

System: Apex

Subsystem: Microsoft Win Pro Embedded 10 LTSC 2016

Subject: Microsoft Win Pro Embedded 10 LTSC 2016 Patch release August 2022

Purpose:

To release the list of validated Microsoft critical patches on the Apex Win Pro Embedded 10 LTSC 2016 systems.

Note 1: If you are using the following versions of APEX, if you are required to update to the latest Hologic validated Windows updates, you must ensure your system(s) have been patched with the "Sentinel Run-Time 7.80.0.0" which can be located in the Windows Control Panel "Program and Features". If needed you can install the patch by ordering "APEX Service Dongle Patch Release" (ASY-12610). Please contact Hologic Support if you have any questions.

System Type	Software Versions
Horizon	APEX v5.6.0.7
Workstation	APEX v5.6.0.7

List of patches/KB's validated for August 2022:

KB5017095 KB5016622 KB5014164
 KB5012170

List of all patches/KB's validated on the Apex System:

KB890830 KB5012596 KB5006669
 KB5017095 KB5012170 KB5005698
 KB5016622 KB5011570 KB5005573
 KB5016058 KB5011495 KB5005043
 KB5015808 KB5010359 KB5004238
 KB5014702 KB5009557 KB5003638
 KB5014164 KB5008873 KB5003197
 KB5014026 KB5008207 KB5001402
 KB5013952 KB5007192 KB5001347

KB5001078	KB4571694	KB4520724
KB5000803	KB4565912	KB4516115
KB5000685	KB4565511	KB4516044
KB5000645	KB4562561	KB4512574
KB4601392	KB4561616	KB4505422
KB4601318	KB4561600	KB4503537
KB4598243	KB4556813	KB4503308
KB4593226	KB4550994	KB4503267
KB4589210	KB4550929	KB4498947
KB4589208	KB4540723	KB4494175
KB4586830	KB4540670	KB4346087
KB4583462	KB4537764	KB4132216
KB4580346	KB4537759	KB4052908
KB4580325	KB4535680	KB4022619
KB4577586	KB4535288	KB4018073
KB4577015	KB4534271	KB2565063
KB4576750	KB4530689	KB2467173