



**Date:** September 15, 2022

**Author:** Software Engineering

**System:** TRIDENT

**Subsystem:** Windows 7 & Software Security Patches

**Subject:** Windows 7 & Software Security Patch release August 2022

---

---

**Purpose:**

To release the list of validated Microsoft critical patches on the TRIDENT System.

**Note:**

- Please check if January 2020 Windows 7 Monthly Patch Disk is installed on system, if not then first Install January 2020 Windows 7 Monthly Patch Disk (SFW-03075).
- Apply SFW-03338 Rev. 004 (Windows 7 ESU License Installation) to the system prior to installation of March 2022 Windows 7 Monthly Patch Disk (SFW-03646) to allow for updates to be installed, due to end of support for Microsoft Windows 7.

**List of patches/KB's validated for August 2022:**

KB5016679

KB5016676

**List of all patches/KB's validated on the TRIDENT System:**

KB976932

KB5014748

KB5012626

KB5016679

KB5014742

KB5012138

KB5016676

KB5014012

KB5012131

KB5016057

KB5013999

KB5011649

KB5015862

KB5013637

KB5011552

KB5015861

KB5012649

KB5011529

|           |           |           |
|-----------|-----------|-----------|
| KB5010422 | KB5000841 | KB4565623 |
| KB5009621 | KB4601363 | KB4565612 |
| KB5008867 | KB4601347 | KB4565539 |
| KB5008859 | KB4600945 | KB4565524 |
| KB5008282 | KB4598289 | KB4565354 |
| KB5007233 | KB4598279 | KB4562030 |
| KB5006749 | KB4597239 | KB4561669 |
| KB5006743 | KB4592510 | KB4561643 |
| KB5006728 | KB4592503 | KB4556843 |
| KB5005633 | KB4592471 | KB4556836 |
| KB5005615 | KB4586827 | KB4556798 |
| KB5005089 | KB4586805 | KB4555449 |
| KB5005088 | KB4585204 | KB4552940 |
| KB5004951 | KB4580970 | KB4552919 |
| KB5004757 | KB4580387 | KB4550965 |
| KB5004378 | KB4580345 | KB4541500 |
| KB5004307 | KB4578963 | KB4538483 |
| KB5004289 | KB4578952 | KB4537829 |
| KB5004120 | KB4577053 | KB4537813 |
| KB5003694 | KB4577051 | KB4536952 |
| KB5003667 | KB4576612 | KB4534314 |
| KB5003233 | KB4575903 | KB4534310 |
| KB5003228 | KB4571729 | KB4533012 |
| KB5001848 | KB4571719 | KB4532945 |
| KB5001392 | KB4570673 | KB4532932 |
| KB5001335 | KB4569775 | KB4531786 |
| KB5000851 | KB4569767 | KB4530734 |

|           |           |           |
|-----------|-----------|-----------|
| KB4530692 | KB4490628 | KB4074837 |
| KB4525235 | KB4489885 | KB4074587 |
| KB4525233 | KB4489878 | KB4057114 |
| KB4523206 | KB4486564 | KB4056897 |
| KB4523206 | KB4486563 | KB4054521 |
| KB4520003 | KB4483458 | KB4048960 |
| KB4519976 | KB4483451 | KB4041678 |
| KB4516655 | KB4483187 | KB4038779 |
| KB4516065 | KB4480970 | KB4034679 |
| KB4516033 | KB4480960 | KB4025337 |
| KB4515854 | KB4480063 | KB4022722 |
| KB4512506 | KB4480055 | KB4022619 |
| KB4512486 | KB4474419 | KB4020322 |
| KB4511516 | KB4471328 | KB4019990 |
| KB4507456 | KB4468323 | KB4019263 |
| KB4507449 | KB4467106 | KB4015546 |
| KB4507004 | KB4462915 | KB4012212 |
| KB4506997 | KB4457145 | KB3212642 |
| KB4505422 | KB4457044 | KB3205394 |
| KB4503292 | KB4457035 | KB3197867 |
| KB4503269 | KB4343899 | KB3192391 |
| KB4499175 | KB4339284 | KB3177467 |
| KB4499164 | KB4338823 | KB3172605 |
| KB4495606 | KB4284867 | KB3161949 |
| KB4495588 | KB4103712 | KB3159398 |
| KB4493472 | KB4096418 | KB3156016 |
| KB4493448 | KB4093108 | KB3150220 |

|           |           |           |
|-----------|-----------|-----------|
| KB3139914 | KB3059317 | KB2943357 |
| KB3139398 | KB3055642 | KB2937610 |
| KB3138910 | KB3046269 | KB2931356 |
| KB3138612 | KB3046017 | KB2912390 |
| KB3133977 | KB3045685 | KB2911501 |
| KB3127220 | KB3042058 | KB2900986 |
| KB3126587 | KB3037574 | KB2894854 |
| KB3124275 | KB3035126 | KB2894844 |
| KB3122648 | KB3030377 | KB2893294 |
| KB3110329 | KB3023215 | KB2888049 |
| KB3109560 | KB3022777 | KB2884256 |
| KB3109103 | KB3021674 | KB2882822 |
| KB3108664 | KB3019978 | KB2871997 |
| KB3108371 | KB3018238 | KB2868038 |
| KB3101722 | KB3011780 | KB2864202 |
| KB3097989 | KB3010788 | KB2862335 |
| KB3092601 | KB2992611 | KB2862330 |
| KB3086255 | KB2991963 | KB2862152 |
| KB3080079 | KB2984972 | KB2861698 |
| KB3078601 | KB2979596 | KB2840631 |
| KB3075220 | KB2978742 | KB2834140 |
| KB3074543 | KB2977292 | KB2813430 |
| KB3072305 | KB2973201 | KB2813347 |
| KB3071756 | KB2973112 | KB2807986 |
| KB3067903 | KB2972211 | KB2789645 |
| KB3061518 | KB2972100 | KB2786081 |
| KB3060716 | KB2968294 | KB2770660 |

|           |           |
|-----------|-----------|
| KB2758857 | KB2560656 |
| KB2742599 | KB2552343 |
| KB2736422 | KB2533623 |
| KB2731771 | KB2533552 |
| KB2729452 | KB2532531 |
| KB2729094 | KB2510531 |
| KB2727528 | KB2506212 |
| KB2719033 | KB2506014 |
| KB2716513 | KB2491683 |
| KB2705219 | KB2479943 |
| KB2698365 | KB2467173 |
| KB2690533 |           |
| KB2685939 |           |
| KB2670838 |           |
| KB2667402 |           |
| KB2656356 |           |
| KB2654428 |           |
| KB2653956 |           |
| KB2639308 |           |
| KB2631813 |           |
| KB2621440 |           |
| KB2620704 |           |
| KB2604115 |           |
| KB2585542 |           |
| KB2579686 |           |
| KB2565063 |           |
| KB2564958 |           |