



Date: September 04, 2025

Author: Software Engineering

System: Path Vision

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release August 2025

Purpose:

To release the list of validated Microsoft critical patches on the Path Vision System.

Note: Core Vision and BioVision are on similar product line as Path Vision and both use the same executable so patches validated on Path Vision system can be installed on Core Visions system and BioVision as well.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.14 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for August 2025:

KB5063877

List of all patches/KB's validated on the Path Vision System:

KB5063877	KB5055519	KB5050008
KB5062557	KB890830	KB5049608
KB5060531	KB5053596	KB5048661
KB5058392	KB5052000	KB5046615

KB5044277	KB5031984	KB5019966
KB5044016	KB5031361	KB5018419
KB5043050	KB5030214	KB5017315
KB5041578	KB5029931	KB5016623
KB5040430	KB5029247	KB5015811
KB5039879	KB5028168	KB5014692
KB5039217	KB5027222	KB5013941
KB5037765	KB5027131	KB5013641
KB5036896	KB5026362	KB5012170
KB5036604	KB5025229	KB5005112
KB5035849	KB5023702	KB4589208
KB5034768	KB5022840	KB4577586
KB5034127	KB5022511	KB4535680
KB5033904	KB5021237	KB4494174
KB5033371	KB5020866	
KB5032196	KB5020627	