



Date: *January 15, 2021*

Author: Software Engineering

System: BREVERA / CORLUMINA

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release
December 2020

Purpose:

To release the list of validated Microsoft critical patches on the BREVERA / CORLUMINA System.

List of patches/KB's validated for December 2020:

KB4592440

List of all patches/KB's validated on the BREVERA / CORLUMINA System:

KB4592440	KB4577667	KB4552924
KB4587735	KB4570333	KB4539571
KB4586793	KB4570332	KB4494174
KB4580325	KB4561600	
KB4577668	KB4557397	