

Date: 01/15/2021

Author: Software Engineering

System: Insight

Subsystem: Microsoft Windows OS 10

Subject: Microsoft Windows 10 OS Patch release December 2020

Purpose:

To release the list of validated Microsoft critical patches on the Insight systems.

System Type	Software Versions
Workstation	Insight v6.1.3

List of patches/KB's validated for December 2020:

KB4593226

List of all patches/KB's validated on the Insight System:

KB4593226	KB4562561	KB4537759
KB4586830	KB4561616	KB4534271
KB4580346	KB4561600	KB4530689
KB4580325	KB4556813	KB4525236
KB4577015	KB4550994	KB4520724
KB4576750	KB4550929	KB4516115
KB4571694	KB4540723	KB4494175
KB4565912	KB4540670	KB4132216
KB4565511	KB4537764	KB2467173