



Date: *January 15, 2021*

Author: Software Engineering

System: UNIFI WORKSPACE SOFTWARE

Subsystem: Windows 10 LTSC & Software Security Patches

Subject: Windows 10 LTSC & Software Security Patch release December 2020

Purpose:

To release the list of validated Microsoft critical patches on the Unifi Workspace Software System.

List of patches/KB's validated for December 2020:

KB4593226	KB4578969	KB4091664
-----------	-----------	-----------

List of all patches/KB's validated on the Unifi Workspace Software System:

KB4593226	KB4562561	KB4494175
KB4586830	KB4561616	KB4346087
KB4580346	KB4561600	KB4132216
KB4580325	KB4556813	KB4091664
KB4578969	KB4550994	KB4052908
KB4577015	KB4550929	KB2565063
KB4576750	KB4549825	KB2467173
KB4565912	KB4540723	
KB4565511	KB4537759	