



Date: January 13, 2026

Author: Software Engineering

System: AWM

Subsystem: Windows Server 2019 & Software Security Patches

Subject: Windows Server 2019 & Software Security Patch release December 2025

Purpose:

To release the list of validated Microsoft critical patches on the AWM System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.21 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

- **.NET Core 2.1 / 2.2**

Hosting Bundle 2.1.30

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for December 2025:

KB5071544

List of all patches/KB's validated on the AWM System:

KB5071544

KB5065225

KB5066143

KB5049608

KB5068791

KB5066586

KB5063877

KB5062557	KB5030214	KB5013952
KB5060531	KB5029931	KB5013641
KB5058392	KB5029376	KB5012596
KB5055519	KB5029247	KB5012170
KB5050533	KB5028168	KB5011570
KB890830	KB5027222	KB5011495
KB5053596	KB5027131	KB5010359
KB5052000	KB5026362	KB5009557
KB5050008	KB5025229	KB5009546
KB5049608	KB5023702	KB5008873
KB5048661	KB5022840	KB5008207
KB5046615	KB5022511	KB5007192
KB5044277	KB5022286	KB5006669
KB5044016	KB5021237	KB5005698
KB5043050	KB5021126	KB5005573
KB5042215	KB5020866	KB5005112
KB5041578	KB5020627	KB5005043
KB5040430	KB5019966	KB5004238
KB5039879	KB5019964	KB5003638
KB5039217	KB5018419	KB5003197
KB5037765	KB5018411	KB5001402
KB5036896	KB5017396	KB5001347
KB5036604	KB5017305	KB5001078
KB5035849	KB5017095	KB5000803
KB5034768	KB5016884	KB5000685
KB5034127	KB5016622	KB4601392
KB5033904	KB5016058	KB4601318
KB5033371	KB5015808	KB4598243
KB5032196	KB5014702	KB4593226
KB5031984	KB5014553	KB4589210
KB5031361	KB5014026	KB4589208

KB4586830

KB4561616

KB4514366

KB4580346

KB4557397

KB4505225

KB4577586

KB4551853

KB4494175

KB4577015

KB4550929

KB4494174

KB4576750

KB4549947

KB4132216

KB4571694

KB4540723

KB2565063

KB4565912

KB4539571

KB2467173

KB4565511

KB4537759

KB4562561

KB4535680