



**Date:** January 13, 2026

**Author:** Software Engineering

**System:** Apex

**Subsystem:** Microsoft Win Pro Embedded 10 LTSP 2016

**Subject:** Microsoft Win Pro Embedded 10 LTSP 2016 Patch release December 2025

---

**Purpose:**

To release the list of validated Microsoft critical patches on the Apex Win Pro Embedded 10 LTSP 2016 systems.

Note: If you are using the following versions of APEX, if you are required to update to the latest Hologic validated Windows updates, you must ensure your system(s) have been patched with the “Sentinel Run-Time 7.80.0.0” which can be located in the Windows Control Panel “Program and Features”. If needed you can install the patch by ordering “APEX Service Dongle Patch Release” (ASY-12610). Please contact Hologic Support if you have any questions.

System Type	Software Versions
Horizon	APEX v5.6.0.7
Workstation	APEX v5.6.0.7

**Latest .NET Versions Validated for the Product:**

**Note:** Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.21 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

**Latest SQLSERVER Versions Validated:**

**Note:** Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2014:**

The latest security updates are included in the Windows Cumulative Update.

**List of patches/KB's validated for December 2025:**

KB5071543

**List of all patches/KB's validated on the Apex System**

KB5071543	KB5065222	
KB5068864	KB5042749	KB5030214
KB5070247	KB5041773	KB5030213
KB5066584	KB5041576	KB5029931
KB5066836	KB5040562	KB5029247
KB5065427	KB5040434	KB5029242
KB5065687	KB5039334	KB5029185
KB5063871	KB5039214	KB5028169
KB5062560	KB5037763	KB5028168
KB5062799	KB5037016	KB5027222
KB5061010	KB5036899	KB5027219
KB5060954	KB5035962	KB5027131
KB5058524	KB5035855	KB5026363
KB5058383	KB5034862	KB5026362
KB5055661	KB5034767	KB5025229
KB5055521	KB5034127	KB5025228
KB5054833	KB5034119	KB5023788
KB890830	KB5033904	KB5023702
KB5054006	KB5033688	KB5023697
KB5053594	KB5033373	KB5022840
KB5052006	KB5033371	KB5022838
KB5050109	KB5032391	KB5022511
KB5049993	KB5032197	KB5022289
KB5048671	KB5032196	KB5022286
KB5046612	KB5031984	KB5021237
KB5044293	KB5031362	KB5021235
KB5043124	KB5031361	KB5021126
KB5043051	KB5030504	KB5021045

KB5020866	KB5010351	KB5001402
KB5020627	KB5009557	KB5001347
KB5019966	KB5009546	KB5001342
KB5019964	KB5008873	KB5001078
KB5018419	KB5008218	KB5000859
KB5018411	KB5008207	KB5000822
KB5017396	KB5007206	KB5000803
KB5017315	KB5007192	KB5000685
KB5017305	KB5006672	KB4601393
KB5017095	KB5006669	KB4601392
KB5016884	KB5005698	KB4601345
KB5016623	KB5005573	KB4601318
KB5016622	KB5005568	KB4601060
KB5016058	KB5005112	KB4598480
KB5015811	KB5005043	KB4598243
KB5015808	KB5005030	KB4598230
KB5014702	KB5004244	KB4593226
KB5014692	KB5004238	KB4592440
KB5014553	KB5003711	KB4589210
KB5014026	KB5003646	KB4589208
KB5013952	KB5003638	KB4587735
KB5013941	KB5003243	KB4586830
KB5013641	KB5003197	KB4586793
KB5012647	KB5003171	KB4583465
KB5012596	KB5002542	KB4580346
KB5012170	KB5002520	KB4580325
KB5012128	KB5002497	KB4577668
KB5011570	KB5002464	KB4577667
KB5011503	KB5002406	KB4577586
KB5011495	KB5001407	KB4577015
KB5010359	KB5001404	KB4576750

KB4571694	KB4539571	KB4512574
KB4570333	KB4538461	KB4512517
KB4570332	KB4537764	KB4511553
KB4566424	KB4537759	KB4509091
KB4565912	KB4535680	KB4507460
KB4565511	KB4534273	KB4506998
KB4565349	KB4534271	KB4505225
KB4562562	KB4532691	KB4503308
KB4562561	KB4530715	KB4498947
KB4561616	KB4530689	KB4494440
KB4561608	KB4525236	KB4494175
KB4561600	KB4523205	KB4494174
KB4558998	KB4523204	KB4493470
KB4558997	KB4521862	KB4489882
KB4557397	KB4521858	KB4487026
KB4556813	KB4520724	KB4485447
KB4552924	KB4519998	KB4346087
KB4551853	KB4519338	KB4132216
KB4550994	KB4516115	KB4091664
KB4550929	KB4516044	KB4022619
KB4549949	KB4514366	KB4018073
KB4549947	KB4512937	KB2565063
KB4540723	KB4512578	KB2467173
KB4540670	KB4512577	