



**Date:** January 13, 2026

**Author:** Software Engineering

**System:** BREVERA / CORLUMINA

**Subsystem:** Windows 10 LTSC 2019 (v1809) & Software Security Patches

**Subject:** Windows 10 LTSC 2019 (v1809) & Software Security Patch release  
November December 2025

**Purpose:**

To release the list of validated Microsoft critical patches on the BREVERA / CORLUMINA System.

**Latest .NET Versions Validated for the Product:**

**Note:** Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.21 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

**Latest SQLSERVER Versions Validated:**

**Note:** Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

**List of patches/KB's validated for December 2025:**

KB5071544

**List of all patches/KB's validated on the BREVERA / CORLUMINA System:**

KB5071544

KB5068791

KB5039879

KB5023788

KB5065225

KB5039217

KB5023697

KB5066143

KB5037765

KB5022838

KB5066586

KB5036896

KB5022511

KB5065428

KB5036604

KB5022289

KB5063877

KB5035849

KB5021235

KB5062557

KB5034768

KB5021126

KB5060531

KB5034127

KB5019964

KB5058392

KB5033904

KB5018411

KB5055519

KB5033371

KB5017396

KB5050533

KB5032196

KB5017305

KB890830

KB5031984

KB5017095

KB5053596

KB5031361

KB5016884

KB5052000

KB5030214

KB5016622

KB5050008

KB5029931

KB5016058

KB5049608

KB5029376

KB5015808

KB5048661

KB5029247

KB5014702

KB5046615

KB5028168

KB5014553

KB5044277

KB5027222

KB5014026

KB5044016

KB5027131

KB5013952

KB5043050

KB5026363

KB5012596

KB5042215

KB5026362

KB5012170

KB5041578

KB5025229

KB5011570

KB5040430

KB5025228

KB5011495

KB5010359	KB4589208	KB4521858
KB5009546	KB4586830	KB4520724
KB5008207	KB4580346	KB4519998
KB5007192	KB4577586	KB4516044
KB5006669	KB4577015	KB4512574
KB5005698	KB4576750	KB4512517
KB5005573	KB4571694	KB4509091
KB5005112	KB4565912	KB4507460
KB5005043	KB4565511	KB4505225
KB5004238	KB4562561	KB4498947
KB5003638	KB4561616	KB4494440
KB5003197	KB4557397	KB4494175
KB5001402	KB4556813	KB4494174
KB5001347	KB4550994	KB4493470
KB5001078	KB4550929	KB4489882
KB5000803	KB4540723	KB4487026
KB5000685	KB4540670	KB4485447
KB4601392	KB4537764	KB4346087
KB4601318	KB4535680	KB4132216
KB4598243	KB4534271	KB4091664
KB4593226	KB4530689	KB2565063
KB4589210	KB4525236	KB2467173