



**Date:** January 13, 2026

**Author:** Software Engineering

**System:** SELENIA DIMENSIONS / 3DIMENSIONS

**Subsystem:** Windows 10 LTSB 2016 (v1607) & Software Security Patches

**Subject:** Windows 10 LTSB 2016 (v1607) & Software Security Patch release  
December 2025

---

**Purpose:**

To release the list of validated Microsoft critical patches on the SELENIA DIMENSIONS / 3DIMENSIONS System.

**Latest .NET Versions Validated for the Product:**

**Note:** Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.21 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

- **.NET Core 2.1 / 2.2**

Hosting Bundle 2.1.30

**Latest SQLSERVER Versions Validated:**

**Note:** Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

**List of patches/KB's validated for December 2025:**

KB5071543

**List of all patches/KB's validated on the SELENIA DIMENSIONS /  
3DIMENSIONS System:**

KB5071543

KB5065225

KB5070247

KB5066836

KB5068864

KB5066584

KB5065427

KB5065687

KB5042215

KB5029242

KB5063871

KB5041773

KB5028169

KB5062560

KB5041576

KB5027219

KB5062799

KB5040562

KB5026363

KB5060954

KB5040434

KB5025228

KB5061010

KB5039334

KB5023788

KB5058524

KB5039214

KB5023697

KB5058383

KB5037763

KB5022838

KB5055661

KB5037016

KB5022289

KB5055521

KB5036899

KB5021235

KB5050533

KB5035962

KB5021126

KB890830

KB5035855

KB5019964

KB5054006

KB5034862

KB5018411

KB5053594

KB5034767

KB5017396

KB5052006

KB5034119

KB5017305

KB5050109

KB5033373

KB5017095

KB5049993

KB5032391

KB5016884

KB5048671

KB5032197

KB5016622

KB5046612

KB5031362

KB5016058

KB5044293

KB5030504

KB5015808

KB5043124

KB5030213

KB5014702

KB5043051

KB5029376

KB5014553

KB5014026	KB4601051	KB4520724
KB5013952	KB4598243	KB4519998
KB5013625	KB4593226	KB4516044
KB5012596	KB4589210	KB4512574
KB5012170	KB4586830	KB4512517
KB5012118	KB4580346	KB4509091
KB5011570	KB4578969	KB4507460
KB5011495	KB4577586	KB4505225
KB5010359	KB4577015	KB4503537
KB5009546	KB4576750	KB4503267
KB5008877	KB4571694	KB4498947
KB5008207	KB4565912	KB4494440
KB5007192	KB4565511	KB4494175
KB5006669	KB4562561	KB4493472
KB5005698	KB4561616	KB4493470
KB5005573	KB4557397	KB4493448
KB5005043	KB4556813	KB4493435
KB5004238	KB4550994	KB4489882
KB5003638	KB4550929	KB4485447
KB5003197	KB4540723	KB4346087
KB5001402	KB4540670	KB4132216
KB5001347	KB4537764	KB4091664
KB5001078	KB4535680	KB2565063
KB5000803	KB4534271	KB2467173
KB5000685	KB4530689	
KB4601392	KB4525236	
KB4601318	KB4521858	