



**Date:** January 13, 2026

**Author:** Software Engineering

**System:** SELENIA DIMENSIONS / 3DIMENSIONS

**Subsystem:** Windows 10 LTSC 2019 (v1809) & Software Security Patches

**Subject:** Windows 10 LTSC 2019 (v1809) & Software Security Patch release  
December 2025

---

**Purpose:**

To release the list of validated Microsoft critical patches on the SELENIA DIMENSIONS / 3DIMENSIONS System.

**Latest .NET Versions Validated for the Product:**

**Note:** Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.21 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

- **.NET Core 2.1 / 2.2**

Hosting Bundle 2.1.30

**Latest SQLSERVER Versions Validated:**

**Note:** Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

**List of patches/KB's validated for December 2025:**

KB5071544

**List of all patches/KB's validated on the SELENIA DIMENSIONS /  
3DIMENSIONS System:**

KB5071544

KB5065225

KB5040430

KB5022840

KB5068791

KB5039879

KB5022504

KB5066137

KB5039217

KB5022286

KB5066586

KB5037765

KB5021237

KB5065428

KB5036896

KB5021126

KB5063877

KB5036604

KB5020874

KB5062557

KB5035849

KB5020615

KB5060531

KB5034768

KB5019966

KB5058392

KB5034127

KB5018419

KB5055519

KB5033904

KB5017315

KB5050533

KB5033371

KB5016884

KB5049615

KB5032196

KB5016623

KB890830

KB5031990

KB5015811

KB5053596

KB5031361

KB5014692

KB5052000

KB5030214

KB5014553

KB5050008

KB5029925

KB5014164

KB5049608

KB5029376

KB5013941

KB5048661

KB5029247

KB5013641

KB5046615

KB5028168

KB5013626

KB5044277

KB5027222

KB5012647

KB5044016

KB5027124

KB5012170

KB5043050

KB5026362

KB5012119

KB5042215

KB5025229

KB5011503

KB5041578

KB5023702

KB5010351

KB5009557	KB4586793	KB2565063
KB5008873	KB4580325	KB2467173
KB5008218	KB4577668	
KB5007206	KB4577667	
KB5006672	KB4577586	
KB5005568	KB4570333	
KB5005112	KB4570332	
KB5005030	KB4566424	
KB5004244	KB4565349	
KB5003711	KB4562562	
KB5003646	KB4561608	
KB5003243	KB4561600	
KB5003171	KB4558998	
KB5001407	KB4558997	
KB5001404	KB4557397	
KB5001342	KB4552924	
KB5000859	KB4551853	
KB5000822	KB4549949	
KB5000685	KB4549947	
KB4601393	KB4539571	
KB4601345	KB4537759	
KB4601060	KB4535680	
KB4598480	KB4514366	
KB4598230	KB4505225	
KB4592440	KB4494174	
KB4589208	KB4022619	
KB4587735	KB4018073	