

Date: 03/09/2021
Author: Software Engineering
System: Apex
Subsystem: Microsoft Windows OS 10
Subject: Microsoft Windows 10 OS Patch release February 2021

Purpose:

To release the list of validated Microsoft critical patches on the Apex systems.

Note 1: If you are using the following versions of APEX, if you are required to update to the latest Hologic validated Windows updates, you must ensure your system(s) have been patched with the “Sentinel Run-Time 7.80.0.0” which can be located in the Windows Control Panel “Program and Features”. If needed you can install the patch by ordering “APEX Service Dongle Patch Release” (ASY-12610). Please contact Hologic Support if you have any questions.

| System Type | Software Versions |
|--------------------|--------------------------|
| Horizon | APEX v5.6.0.7 |
| Workstation | APEX v5.6.0.7 |

List of patches/KB’s validated for February 2021:

KB4601392 KB4601318

List of all patches/KB’s validated on the Apex System:

| | | |
|-----------|-----------|-----------|
| KB4601392 | KB4571694 | KB4540723 |
| KB4601318 | KB4565912 | KB4540670 |
| KB4598243 | KB4565511 | KB4537764 |
| KB4593226 | KB4562561 | KB4537759 |
| KB4586830 | KB4561616 | KB4535680 |
| KB4580346 | KB4561600 | KB4535288 |
| KB4580325 | KB4556813 | KB4534271 |
| KB4577015 | KB4550994 | KB4530689 |
| KB4576750 | KB4550929 | KB4520724 |

KB4516115

KB4503308

KB4132216

KB4516044

KB4503267

KB4022619

KB4512574

KB4498947

KB4018073

KB4505422

KB4494175

KB2565063

KB4503537

KB4346087

KB2467173