



**Date:** *March 09, 2021*

**Author:** Software Engineering

**System:** BREVERA / CORLUMINA

**Subsystem:** Windows 10 LTSC 2019 (v1809) & Software Security Patches

**Subject:** Windows 10 LTSC 2019 (v1809) & Software Security Patch release  
February 2021

---

**Purpose:**

To release the list of validated Microsoft critical patches on the BREVERA / CORLUMINA System.

**List of patches/KB's validated for February 2021:**

KB4601393                      KB4601345

**List of all patches/KB's validated on the BREVERA / CORLUMINA System:**

KB4601393	KB4586793	KB4561600
KB4601345	KB4580325	KB4557397
KB4598480	KB4577668	KB4552924
KB4598230	KB4577667	KB4539571
KB4592440	KB4570333	KB4535680
KB4587735	KB4570332	KB4494174