

Date: 03/04/2022

Author: Software Engineering

System: Apex

Subsystem: Microsoft Win Pro Embedded 10 LTSB 2016

Subject: Microsoft Win Pro Embedded 10 LTSB 2016 Patch release February 2022

To release the list of validated Microsoft critical patches on the Apex Win Pro Embedded 10 LTSB 2016 systems.

Note 1: If you are using the following versions of APEX, if you are required to update to the latest Hologic validated Windows updates, you must ensure your system(s) have been patched with the “Sentinel Run-Time 7.80.0.0” which can be located in the Windows Control Panel “Program and Features”. If needed you can install the patch by ordering “APEX Service Dongle Patch Release” (ASY-12610). Please contact Hologic Support if you have any questions.

System Type	Software Versions
Horizon	APEX v5.6.0.7
Workstation	APEX v5.6.0.7

List of patches/KB’s validated for February 2022:

KB5010359 KB5000685 KB4589210

List of all patches/KB’s validated on the Apex System:

KB890830 KB5005043 KB5000645
KB5010359 KB5004238 KB4601392
KB5009557 KB5003638 KB4601318
KB5008873 KB5003197 KB4598243
KB5008207 KB5001402 KB4593226
KB5007192 KB5001347 KB4589210
KB5006669 KB5001078 KB4589208
KB5005698 KB5000803 KB4586830
KB5005573 KB5000685 KB4583462

KB4580346	KB4494175
KB4580325	KB4346087
KB4577586	KB4132216
KB4577015	KB4052908
KB4576750	KB4022619
KB4571694	KB4018073
KB4565912	KB2565063
KB4565511	KB2467173
KB4562561	
KB4561616	
KB4561600	
KB4556813	
KB4550994	
KB4550929	
KB4540723	
KB4540670	
KB4537764	
KB4537759	
KB4535680	
KB4535288	
KB4534271	
KB4530689	
KB4520724	
KB4516115	
KB4516044	
KB4512574	
KB4505422	
KB4503537	
KB4503308	
KB4503267	
KB4498947	