



**Date:** March 12, 2026

**Author:** Software Engineering

**System:** Apex

**Subsystem:** Microsoft Win Pro Embedded 10 LTSP 2016

**Subject:** Microsoft Win Pro Embedded 10 LTSP 2016 Patch release February 2026

---

**Purpose:**

To release the list of validated Microsoft critical patches on the Apex Win Pro Embedded 10 LTSP 2016 systems.

Note: If you are using the following versions of APEX, if you are required to update to the latest Hologic validated Windows updates, you must ensure your system(s) have been patched with the “Sentinel Run-Time 7.80.0.0” which can be located in the Windows Control Panel “Program and Features”. If needed you can install the patch by ordering “APEX Service Dongle Patch Release” (ASY-12610). Please contact Hologic Support if you have any questions.

System Type	Software Versions
Horizon	APEX v5.6.0.7
Workstation	APEX v5.6.0.7

**Latest .NET Versions Validated for the Product:**

**Note:** Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.21 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

**Latest SQLSERVER Versions Validated:**

**Note:** Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2014:**

The latest security updates are included in the Windows Cumulative Update.

**List of patches/KB's validated for February 2026:**

KB5075999

KB5075902

**List of all patches/KB's validated on the Apex System**

KB5075902

KB5073447

KB5071543

KB5075999

KB5073722

KB5065222

KB5068864

KB5043124

KB5031984

KB5070247

KB5043051

KB5031362

KB5066584

KB5042749

KB5031361

KB5066836

KB5041773

KB5030504

KB5065427

KB5041576

KB5030214

KB5065687

KB5040562

KB5030213

KB5063871

KB5040434

KB5029931

KB5062560

KB5039334

KB5029247

KB5062799

KB5039214

KB5029242

KB5061010

KB5037763

KB5029185

KB5060954

KB5037016

KB5028169

KB5058524

KB5036899

KB5028168

KB5058383

KB5035962

KB5027222

KB5055661

KB5035855

KB5027219

KB5055521

KB5034862

KB5027131

KB5054833

KB5034767

KB5026363

KB890830

KB5034127

KB5026362

KB5054006

KB5034119

KB5025229

KB5053594

KB5033904

KB5025228

KB5052006

KB5033688

KB5023788

KB5050109

KB5033373

KB5023702

KB5049993

KB5033371

KB5023697

KB5048671

KB5032391

KB5022840

KB5046612

KB5032197

KB5022838

KB5044293

KB5032196

KB5022511

KB5022289	KB5012170	KB5002520
KB5022286	KB5012128	KB5002497
KB5021237	KB5011570	KB5002464
KB5021235	KB5011503	KB5002406
KB5021126	KB5011495	KB5001407
KB5021045	KB5010359	KB5001404
KB5020866	KB5010351	KB5001402
KB5020627	KB5009557	KB5001347
KB5019966	KB5009546	KB5001342
KB5019964	KB5008873	KB5001078
KB5018419	KB5008218	KB5000859
KB5018411	KB5008207	KB5000822
KB5017396	KB5007206	KB5000803
KB5017315	KB5007192	KB5000685
KB5017305	KB5006672	KB4601393
KB5017095	KB5006669	KB4601392
KB5016884	KB5005698	KB4601345
KB5016623	KB5005573	KB4601318
KB5016622	KB5005568	KB4601060
KB5016058	KB5005112	KB4598480
KB5015811	KB5005043	KB4598243
KB5015808	KB5005030	KB4598230
KB5014702	KB5004244	KB4593226
KB5014692	KB5004238	KB4592440
KB5014553	KB5003711	KB4589210
KB5014026	KB5003646	KB4589208
KB5013952	KB5003638	KB4587735
KB5013941	KB5003243	KB4586830
KB5013641	KB5003197	KB4586793
KB5012647	KB5003171	KB4583465
KB5012596	KB5002542	KB4580346

KB4580325	KB4549949	KB4512578
KB4577668	KB4549947	KB4512577
KB4577667	KB4540723	KB4512574
KB4577586	KB4540670	KB4512517
KB4577015	KB4539571	KB4511553
KB4576750	KB4538461	KB4509091
KB4571694	KB4537764	KB4507460
KB4570333	KB4537759	KB4506998
KB4570332	KB4535680	KB4505225
KB4566424	KB4534273	KB4503308
KB4565912	KB4534271	KB4498947
KB4565511	KB4532691	KB4494440
KB4565349	KB4530715	KB4494175
KB4562562	KB4530689	KB4494174
KB4562561	KB4525236	KB4493470
KB4561616	KB4523205	KB4489882
KB4561608	KB4523204	KB4487026
KB4561600	KB4521862	KB4485447
KB4558998	KB4521858	KB4346087
KB4558997	KB4520724	KB4132216
KB4557397	KB4519998	KB4091664
KB4556813	KB4519338	KB4022619
KB4552924	KB4516115	KB4018073
KB4551853	KB4516044	KB2565063
KB4550994	KB4514366	KB2467173
KB4550929	KB4512937	