



Date: March 13, 2026

Author: Software Engineering

System: BREVERA / CORLUMINA

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release
February 2026

Purpose:

To release the list of validated Microsoft critical patches on the BREVERA / CORLUMINA System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.21 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for February 2026:

KB5075904

List of all patches/KB's validated on the BREVERA / CORLUMINA System:

KB5075904	KB5073723	KB5071544
KB5068791	KB5040430	KB5025229
KB5065225	KB5039879	KB5025228
KB5066143	KB5039217	KB5023788
KB5066586	KB5037765	KB5023697
KB5065428	KB5036896	KB5022838
KB5063877	KB5036604	KB5022511
KB5062557	KB5035849	KB5022289
KB5060531	KB5034768	KB5021235
KB5058392	KB5034127	KB5021126
KB5055519	KB5033904	KB5019964
KB5050533	KB5033371	KB5018411
KB890830	KB5032196	KB5017396
KB5053596	KB5031984	KB5017305
KB5052000	KB5031361	KB5017095
KB5050008	KB5030214	KB5016884
KB5049608	KB5029931	KB5016622
KB5048661	KB5029376	KB5016058
KB5046615	KB5029247	KB5015808
KB5044277	KB5028168	KB5014702
KB5044016	KB5027222	KB5014553
KB5043050	KB5027131	KB5014026
KB5042215	KB5026363	KB5013952
KB5041578	KB5026362	KB5012596

KB5012170	KB4593226	KB4525236
KB5011570	KB4589210	KB4521858
KB5011495	KB4589208	KB4520724
KB5010359	KB4586830	KB4519998
KB5009546	KB4580346	KB4516044
KB5008207	KB4577586	KB4512574
KB5007192	KB4577015	KB4512517
KB5006669	KB4576750	KB4509091
KB5005698	KB4571694	KB4507460
KB5005573	KB4565912	KB4505225
KB5005112	KB4565511	KB4498947
KB5005043	KB4562561	KB4494440
KB5004238	KB4561616	KB4494175
KB5003638	KB4557397	KB4494174
KB5003197	KB4556813	KB4493470
KB5001402	KB4550994	KB4489882
KB5001347	KB4550929	KB4487026
KB5001078	KB4540723	KB4485447
KB5000803	KB4540670	KB4346087
KB5000685	KB4537764	KB4132216
KB4601392	KB4535680	KB4091664
KB4601318	KB4534271	KB2565063
KB4598243	KB4530689	KB2467173