



Date: March 13, 2026

Author: Software Engineering

System: Faxitron Core Vision

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release February 2026

Purpose:

To release the list of validated Microsoft critical patches on the Faxitron Core Vision System.

Note: Path Vision and BioVision are on similar product line as Core Vision, and both use the same executable so patches validated on Core Vision system can be installed on Path Visions system and BioVision as well.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.21 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for February 2026:

KB5075904

List of all patches/KB's validated on the Faxitron Core Vision System:

KB5075904	KB5066143	KB5062557
KB5073723	KB5066586	KB5060531
KB5071544	KB5065428	KB5058392
KB5068791	KB5063877	KB5055519

KB890830	KB5034768	KB5021237
KB5053596	KB5034127	KB5020866
KB5052000	KB5033904	KB5020627
KB5050008	KB5033371	KB5019966
KB5049608	KB5032196	KB5018419
KB5048661	KB5031984	KB5017315
KB5046615	KB5031361	KB5016623
KB5044277	KB5030214	KB5015811
KB5044016	KB5029931	KB5014692
KB5043050	KB5029247	KB5013941
KB5041578	KB5028168	KB5013641
KB5040430	KB5027222	KB5012170
KB5039879	KB5027131	KB5005112
KB5039217	KB5026362	KB4589208
KB5037765	KB5025229	KB4577586
KB5036896	KB5023702	KB4535680
KB5036604	KB5022840	KB4494174
KB5035849	KB5022511	