



Date: March 13, 2026

Author: Software Engineering

System: Faxitron Path Vision

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release February 2026

Purpose:

To release the list of validated Microsoft critical patches on the Faxitron Path Vision System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.21 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for February 2026:

KB5075904

List of all patches/KB's validated on the Faxitron Path Vision System:

KB5075904	KB5065428	KB890830
KB5073723	KB5063877	KB5053596
KB5071544	KB5062557	KB5052000
KB5068791	KB5060531	KB5050008
KB5066143	KB5058392	KB5049608
KB5066586	KB5055519	KB5048661

KB5046615	KB5032196	KB5020627
KB5044277	KB5031984	KB5019966
KB5044016	KB5031361	KB5018419
KB5043050	KB5030214	KB5017315
KB5041578	KB5029931	KB5016623
KB5040430	KB5029247	KB5015811
KB5039879	KB5028168	KB5014692
KB5039217	KB5027222	KB5013941
KB5037765	KB5027131	KB5013641
KB5036896	KB5026362	KB5012170
KB5036604	KB5025229	KB5005112
KB5035849	KB5023702	KB4589208
KB5034768	KB5022840	KB4577586
KB5034127	KB5022511	KB4535680
KB5033904	KB5021237	KB4494174
KB5033371	KB5020866	