

**Date:** March 13, 2026

**Author:** Software Engineering

**System:** InSight

**Subsystem:** Microsoft Windows 10 Enterprise 2019 LTSC

**Subject:** Microsoft Windows 10 Enterprise 2019 LTSC Patch release February 2026

---

**Purpose:**

To release the list of validated Microsoft critical patches on the Insights systems.

System Type	Software Versions
Workstation	InSight v6.1.4

**Latest .NET Versions Validated for the Product:**

**Note:** Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.21 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

**List of patches/KB's validated for February 2026:**

KB5075904

**List of all patches/KB's validated on the Insight System:**

KB5075904	KB5058392	KB5044016
KB5073723	KB5055519	KB5043050
KB5071544	KB890830	KB5041578
KB5068791	KB5053596	KB5040430
KB5066143	KB5052000	KB5039879
KB5066586	KB5050008	KB5039217
KB5065428	KB5049608	KB5037765
KB5063877	KB5048661	KB5036896
KB5062557	KB5046615	KB5036604
KB5060531	KB5044277	KB5035849

KB5034768  
KB5034127  
KB5033904  
KB5033371  
KB5032196  
KB5031984  
KB5031361  
KB5030214  
KB5029931  
KB5029247  
KB5028168  
KB5027222  
KB5027131  
KB5026362  
KB5025229  
KB5023702  
KB5022840  
KB5022511  
KB5012170  
KB5005112  
KB5004424  
KB4589208  
KB4577586