



Date: March 13, 2026

Author: Software Engineering

System: SELENIA DIMENSIONS / 3DIMENSIONS

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release
February 2026

Purpose:

To release the list of validated Microsoft critical patches on the SELENIA DIMENSIONS / 3DIMENSIONS System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.21 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

- **.NET Core 2.1 / 2.2**

Hosting Bundle 2.1.30

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for February 2026:

KB5075904

**List of all patches/KB's validated on the SELENIA DIMENSIONS /
3DIMENSIONS System:**

KB5075904	KB5044277	KB5028168
KB5066143	KB5044016	KB5027222
KB5073723	KB5043050	KB5027124
KB5071544	KB5042215	KB5026362
KB5065225	KB5041578	KB5025229
KB5068791	KB5040430	KB5023702
KB5066137	KB5039879	KB5022840
KB5066586	KB5039217	KB5022504
KB5065428	KB5037765	KB5022286
KB5063877	KB5036896	KB5021237
KB5062557	KB5036604	KB5021126
KB5060531	KB5035849	KB5020874
KB5058392	KB5034768	KB5020615
KB5055519	KB5034127	KB5019966
KB5050533	KB5033904	KB5018419
KB5049615	KB5033371	KB5017315
KB890830	KB5032196	KB5016884
KB5053596	KB5031990	KB5016623
KB5052000	KB5031361	KB5015811
KB5050008	KB5030214	KB5014692
KB5049608	KB5029925	KB5014553
KB5048661	KB5029376	KB5014164
KB5046615	KB5029247	KB5013941

KB5013641	KB5001404	KB4565349
KB5013626	KB5001342	KB4562562
KB5012647	KB5000859	KB4561608
KB5012170	KB5000822	KB4561600
KB5012119	KB5000685	KB4558998
KB5011503	KB4601393	KB4558997
KB5010351	KB4601345	KB4557397
KB5009557	KB4601060	KB4552924
KB5008873	KB4598480	KB4551853
KB5008218	KB4598230	KB4549949
KB5007206	KB4592440	KB4549947
KB5006672	KB4589208	KB4539571
KB5005568	KB4587735	KB4537759
KB5005112	KB4586793	KB4535680
KB5005030	KB4580325	KB4514366
KB5004244	KB4577668	KB4505225
KB5003711	KB4577667	KB4494174
KB5003646	KB4577586	KB4022619
KB5003243	KB4570333	KB4018073
KB5003171	KB4570332	KB2565063
KB5001407	KB4566424	KB2467173