



Date: March 13, 2026

Author: Software Engineering

System: SecurView Manager – Windows Server 2016

Subsystem: Windows 10 Server 2016 & Software Security Patches

Subject: Windows 10 Server 2016 & Software Security Patch release February 2026

Purpose:

To release the list of validated Microsoft critical patches on the SecurView Manager – Windows Server 2016 System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.21 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for February 2026:

KB5075999

KB5075902

List of all patches/KB's validated on the SecurView Manager – Windows Server 2016 System:

KB5075902

KB5068864

KB5065687

KB5075999

KB5070247

KB5063871

KB5073722

KB5066584

KB5062560

KB5073447

KB5066836

KB5062799

KB5071543

KB5065427

KB5061010

KB5060954	KB5034767
KB5058524	KB5034119
KB5055661	KB5033373
KB5055521	KB5032391
KB890830	KB5032197
KB5054006	KB5031362
KB5053594	KB5030504
KB5052006	KB5030213
KB5050109	KB5029242
KB5049993	KB5028169
KB5048671	KB5027219
KB5046612	KB5026363
KB5044293	KB5025228
KB5043124	KB5023788
KB5043051	KB5023697
KB5041773	KB5022838
KB5041576	KB5022289
KB5040562	KB5017396
KB5040434	KB5012170
KB5039334	KB5005698
KB5039214	KB4589210
KB5037763	KB4494175
KB5037016	KB4132216
KB5036899	
KB5035962	
KB5035855	
KB5034862	