



Date: 01/20/2020

Author: Software Engineering

System: Apex

Subsystem: Microsoft Windows OS 7

Subject: Microsoft Windows 7 OS Patch release January 2020

Purpose:

To release the list of validated Microsoft critical patches on the Apex systems.

Note: If you are using the following versions of APEX, if you are required to update to the latest Hologic validated Windows updates, you must ensure your system(s) have been patched with the “Sentinel Run-Time 7.80.0.0” which can be located in the Windows Control Panel “Program and Features”. If needed you can install the patch by ordering “APEX Service Dongle Patch Release” (ASY-12610). Please contact Hologic Support if you have any questions.

System Type	Software Versions
Horizon	APEX v5.6.0.7
Workstation	APEX v5.6.0.7

List of patches/KB’s validated for January 2020:

KB4536952 KB4534251 KB4022619
KB4534314 KB4532945
KB4534310 KB4532932

List of all patches/KB’s validated on the Apex System:

KB976932 KB4532945 KB4519108
KB4536952 KB4532932 KB4516033
KB4534314 KB4530692 KB4512486
KB4534310 KB4525233 KB4507704
KB4534251 KB4520003 KB4507456

KB4507004	KB4038779	KB3109103
KB4503269	KB4034679	KB3108664
KB4501226	KB4025337	KB3108371
KB4499175	KB4022722	KB3101722
KB4493448	KB4022619	KB3097989
KB4490628	KB4019990	KB3093513
KB4490128	KB4019263	KB3092601
KB4489885	KB4018073	KB3086255
KB4486564	KB4015546	KB3080079
KB4486459	KB4012212	KB3075226
KB4483187	KB3205394	KB3075220
KB4480960	KB3197867	KB3074543
KB4474419	KB3192391	KB3072305
KB4471328	KB3172605	KB3071756
KB4467106	KB3161949	KB3067903
KB4462915	KB3159398	KB3061518
KB4457145	KB3156016	KB3060716
KB4457044	KB3155178	KB3059317
KB4457035	KB3150220	KB3046269
KB4343899	KB3139914	KB3046017
KB4338823	KB3139398	KB3045685
KB4284867	KB3138910	KB3042058
KB4103712	KB3138612	KB3037574
KB4096418	KB3133977	KB3035126
KB4093108	KB3127220	KB3031432
KB4074587	KB3126587	KB3030377
KB4057116	KB3124275	KB3023215
KB4056897	KB3122648	KB3022777
KB4054521	KB3115858	KB3021674
KB4048960	KB3110329	KB3020388
KB4041678	KB3109560	KB3019978

KB3011780	KB2862152	KB2620704
KB3010788	KB2861698	KB2604115
KB3004375	KB2857650	KB2585542
KB2992611	KB2852386	KB2579686
KB2991963	KB2847927	KB2574819
KB2984976	KB2840631	KB2565063
KB2984972	KB2834140	KB2564958
KB2978742	KB2830477	KB2560656
KB2977292	KB2813430	KB2552343
KB2973201	KB2813347	KB2533623
KB2973112	KB2807986	KB2533552
KB2972211	KB2786081	KB2532531
KB2972100	KB2770660	KB2506212
KB2968294	KB2758857	KB2491683
KB2943357	KB2742599	KB2479943
KB2937610	KB2736422	KB2467173
KB2931356	KB2731771	
KB2923545	KB2729094	
KB2912390	KB2727528	
KB2911501	KB2705219	
KB2900986	KB2698365	
KB2894854	KB2690533	
KB2894844	KB2685939	
KB2893294	KB2670838	
KB2888049	KB2667402	
KB2884256	KB2656356	
KB2882822	KB2654428	
KB2871997	KB2653956	
KB2868038	KB2639308	
KB2864202	KB2631813	
KB2862330	KB2621440	