

**Date:** *January 27, 2021*

**Author:** Software Engineering

**System:** AWM

**Subsystem:** Windows Server 2016 & Software Security Patches

**Subject:** Windows Server 2016 & Software Security Patch release January 2021

---

**Purpose:**

To release the list of validated Microsoft critical patches on the AWM System.

**List of patches/KB's validated for January 2021:**

KB4598243                      KB4535680

**List of all patches/KB's validated on the AWM System:**

KB4598243                      KB4562561                      KB4535680

KB4593226                      KB4561616                      KB4514366

KB4586830                      KB4557397                      KB4505225

KB4580346                      KB4551853                      KB4494175

KB4577015                      KB4550929                      KB4494174

KB4576750                      KB4549947                      KB4132216

KB4571694                      KB4540723                      KB2565063

KB4565912                      KB4539571                      KB2467173

KB4565511                      KB4537759