

Date: 01/31/2022

Author: Software Engineering

System: Apex

Subsystem: Microsoft Windows OS 7

Subject: Microsoft Windows 7 OS Patch release January 2022

To release the list of validated Microsoft critical patches on the Apex systems.

Note 1: If you are using the following versions of APEX, if you are required to update to the latest Hologic validated Windows updates, you must ensure your system(s) have been patched with the “Sentinel Run-Time 7.80.0.0” which can be located in the Windows Control Panel “Program and Features”. If needed you can install the patch by ordering “APEX Service Dongle Patch Release” (ASY-12610). Please contact Hologic Support if you have any questions.

System Type	Software Versions
Horizon	APEX v5.6.0.7
Workstation	APEX v5.6.0.7

Note 2:

- Please check if January 2020 Windows 7 Monthly Patch Disk is installed on system, if not then first Install January 2020 Windows 7 Monthly Patch Disk (SFW-03075).
- Apply SFW-03338 Rev. 001 (Windows 7 ESU License Installation) to the system prior to installation of January 2021 Windows 7 Monthly Patch Disk (SFW-03375) to allow for updates to be installed, due to end of support for Microsoft Windows 7.
- Apply SFW-03338 Rev. 003 (Windows 7 ESU License Installation) to the system prior to installation of April 2021 Windows 7 Monthly Patch Disk (SFW-03430) to allow for updates to be installed, due to end of support for Microsoft Windows 7.

List of patches/KB's validated for January 2022:

KB5009621 KB5008867

KB5009610 KB5008859

List of all patches/KB's validated on the Apex System:

KB976932 KB5009621 KB5009610

KB5008867	KB5001335	KB4537813
KB5008859	KB5000851	KB4536952
KB5008282	KB5000800	KB4534314
KB5008244	KB4601363	KB4534310
KB5007236	KB4601275	KB4534251
KB5007233	KB4600945	KB4532945
KB5006749	KB4598289	KB4532932
KB5006743	KB4598279	KB4530692
KB5006728	KB4597239	KB4525233
KB5006671	KB4592510	KB4520003
KB5005633	KB4592503	KB4519108
KB5005615	KB4592471	KB4516033
KB5005563	KB4586827	KB4512486
KB5005089	KB4586805	KB4507704
KB5005088	KB4586768	KB4507456
KB5005036	KB4585204	KB4507004
KB5004951	KB4580970	KB4503269
KB5004757	KB4580387	KB4501226
KB5004378	KB4578952	KB4499175
KB5004307	KB4578623	KB4493448
KB5004289	KB4577053	KB4490628
KB5004233	KB4575903	KB4490128
KB5004120	KB4571719	KB4489885
KB5003694	KB4566371	KB4486564
KB5003667	KB4565539	KB4486459
KB5003636	KB4561669	KB4483187
KB5003233	KB4557900	KB4480960
KB5003228	KB4556843	KB4474419
KB5003165	KB4550965	KB4471328
KB5001848	KB4541500	KB4467106
KB5001392	KB4537829	KB4462915

KB4457145	KB3156016	KB3060716
KB4457044	KB3155178	KB3059317
KB4457035	KB3150220	KB3046269
KB4343899	KB3139914	KB3046017
KB4338823	KB3139398	KB3045685
KB4284867	KB3138910	KB3042058
KB4103712	KB3138612	KB3037574
KB4096418	KB3133977	KB3035126
KB4093108	KB3127220	KB3031432
KB4074587	KB3126587	KB3030377
KB4057116	KB3124275	KB3023215
KB4056897	KB3122648	KB3022777
KB4054521	KB3115858	KB3021674
KB4048960	KB3110329	KB3020388
KB4041678	KB3109560	KB3019978
KB4038779	KB3109103	KB3011780
KB4034679	KB3108664	KB3010788
KB4025337	KB3108371	KB3004375
KB4022722	KB3101722	KB2992611
KB4022619	KB3097989	KB2991963
KB4019990	KB3093513	KB2984976
KB4019263	KB3092601	KB2984972
KB4018073	KB3086255	KB2978742
KB4015546	KB3080079	KB2977292
KB4012212	KB3075226	KB2973201
KB3205394	KB3075220	KB2973112
KB3197867	KB3074543	KB2972211
KB3192391	KB3072305	KB2972100
KB3172605	KB3071756	KB2968294
KB3161949	KB3067903	KB2943357
KB3159398	KB3061518	KB2937610

KB2931356	KB2731771	KB2467173
KB2923545	KB2729094	
KB2912390	KB2727528	
KB2911501	KB2705219	
KB2900986	KB2698365	
KB2894854	KB2690533	
KB2894844	KB2685939	
KB2893294	KB2670838	
KB2888049	KB2667402	
KB2884256	KB2656356	
KB2882822	KB2654428	
KB2871997	KB2653956	
KB2868038	KB2639308	
KB2864202	KB2631813	
KB2862330	KB2621440	
KB2862152	KB2620704	
KB2861698	KB2604115	
KB2857650	KB2585542	
KB2852386	KB2579686	
KB2847927	KB2574819	
KB2840631	KB2565063	
KB2834140	KB2564958	
KB2830477	KB2560656	
KB2813430	KB2552343	
KB2813347	KB2533623	
KB2807986	KB2533552	
KB2786081	KB2532531	
KB2770660	KB2506212	
KB2758857	KB2491683	
KB2742599	KB2479943	
KB2736422	KB2467173	