



Date: February 8, 2023

Author: Software Engineering

System: Apex

Subsystem: Microsoft Win Pro Embedded 10 LTSP 2016

Subject: Microsoft Win Pro Embedded 10 LTSP 2016 Patch release January 2023

Purpose:

To release the list of validated Microsoft critical patches on the Apex Win Pro Embedded 10 LTSP 2016 systems.

Note 1: If you are using the following versions of APEX, if you are required to update to the latest Hologic validated Windows updates, you must ensure your system(s) have been patched with the "Sentinel Run-Time 7.80.0.0" which can be located in the Windows Control Panel "Program and Features". If needed you can install the patch by ordering "APEX Service Dongle Patch Release" (ASY-12610). Please contact Hologic Support if you have any questions.

System Type	Software Versions
Horizon	APEX v5.6.0.7
Workstation	APEX v5.6.0.7

List of patches/KB's validated for January 2023:

KB5022289 KB4583465

List of all patches/KB's validated on the Apex System:

KB890830 KB5016884 KB5011570
KB5022289 KB5016622 KB5011495
KB5021235 KB5016058 KB5010359
KB5020627 KB5015808 KB5009557
KB5019966 KB5014702 KB5008873
KB5019964 KB5014164 KB5008207
KB5018411 KB5014026 KB5007192
KB5017396 KB5013952 KB5006669
KB5017305 KB5012596 KB5005698
KB5017095 KB5012170 KB5005573

KB5005112	KB4561600
KB5005043	KB4556813
KB5004238	KB4550994
KB5003638	KB4550929
KB5003197	KB4540723
KB5002223	KB4540670
KB5001402	KB4537764
KB5001347	KB4537759
KB5001078	KB4535680
KB5000803	KB4535288
KB5000685	KB4534271
KB5000645	KB4530689
KB4601392	KB4520724
KB4601318	KB4516115
KB4598243	KB4516044
KB4593226	KB4512574
KB4589210	KB4505422
KB4589208	KB4503537
KB4586830	KB4503308
KB4583465	KB4503267
KB4583462	KB4498947
KB4580346	KB4494175
KB4580325	KB4346087
KB4577586	KB4132216
KB4577015	KB4052908
KB4576750	KB4022619
KB4571694	KB4018073
KB4565912	KB2565063
KB4565511	KB2467173
KB4562561	
KB4561616	