



The Science of Sure

**Date:** February 8, 2023

**Author:** Software Engineering

**System:** Core Vision

**Subsystem:** Windows 10 LTSC 2019 (v1809) & Software Security Patches

**Subject:** Windows 10 LTSC 2019 (v1809) & Software Security Patch release January 2023

---

**Purpose:**

To release the list of validated Microsoft critical patches on the Core Vision System.

**List of patches/KB's validated for January 2023:**

Hologic is investigating issues on Faxitron Core Vision that may be related to installation of the latest January 2023 security patches. These patches are currently not validated or recommended for installation. Latest patches tested and recommended for installation are December 2022 patches. The corresponding monthly patch disk for January 2023, SFW-03800, is not recommended for use on Core Vision product systems at this time.

**List of all patches/KB's validated on the Core Vision System:**

KB5021237	KB5016623	KB5005112
KB5020866	KB5015811	KB4589208
KB5020627	KB5014692	KB4577586
KB5019966	KB5013941	KB4535680
KB5018419	KB5013641	
KB5017315	KB5012170	