



**Date:** February 12, 2026

**Author:** Software Engineering

**System:** AWM

**Subsystem:** Windows Server 2016 & Software Security Patches

**Subject:** Windows Server 2016 & Software Security Patch release January 2026

---

**Purpose:**

To release the list of validated Microsoft critical patches on the AWM System.

**Latest .NET Versions Validated for the Product:**

**Note:** Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.21 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

- **.NET Core 2.1 / 2.2**

Hosting Bundle 2.1.30

**Latest SQLSERVER Versions Validated:**

**Note:** Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

**List of patches/KB's validated for January 2026:**

KB5073722

KB5073447

**List of all patches/KB's validated on the AWM System:**

KB5073447

KB5073722

KB5071543

KB5065225

KB5070247

KB5066836

KB5068864

KB5066584

KB5065687

KB5065427	KB5037763	KB5016884
KB5063871	KB5037016	KB5016622
KB5062560	KB5036899	KB5016058
KB5062799	KB5035962	KB5015808
KB5060954	KB5035855	KB5014702
KB5061010	KB5034862	KB5014553
KB5058392	KB5034767	KB5014026
KB5049608	KB5034119	KB5013952
KB5055661	KB5033373	KB5012596
KB5055521	KB5032391	KB5012170
KB5050533	KB5032197	KB5011570
KB890830	KB5031362	KB5011495
KB5054006	KB5030504	KB5010359
KB5053594	KB5030213	KB5009557
KB5052006	KB5029376	KB5009546
KB5050109	KB5029242	KB5008873
KB5049993	KB5028169	KB5008207
KB5048671	KB5027219	KB5007192
KB5046612	KB5026363	KB5006669
KB5044293	KB5025228	KB5005698
KB5043124	KB5023788	KB5005573
KB5043051	KB5023697	KB5005112
KB5042215	KB5022838	KB5005043
KB5041773	KB5022289	KB5004238
KB5041578	KB5021235	KB5003638
KB5041576	KB5021126	KB5003197
KB5040562	KB5019964	KB5001402
KB5040434	KB5018411	KB5001347
KB5039879	KB5017396	KB5001078
KB5039334	KB5017305	KB5000803
KB5039214	KB5017095	KB5000685

KB4601392	KB4571694	KB4537759
KB4601318	KB4565912	KB4535680
KB4598243	KB4565511	KB4514366
KB4593226	KB4562561	KB4505225
KB4589210	KB4561616	KB4494175
KB4589208	KB4557397	KB4494174
KB4586830	KB4551853	KB4132216
KB4580346	KB4550929	KB4132214
KB4577586	KB4549947	KB2565063
KB4577015	KB4540723	KB2467173
KB4576750	KB4539571	