



**Date:** February 10, 2026

**Author:** Software Engineering

**System:** Affirm Prone

**Subsystem:** Windows 10 & Software Security Patches

**Subject:** Windows 10 & Software Security Patch release January 2026

**Purpose:**

To release the list of validated Microsoft critical patches on the Affirm Prone System.

**Latest .NET Versions Validated for the Product:**

**Note:** Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.21 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

**Latest SQLSERVER Versions Validated:**

**Note:** Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

**List of patches/KB's validated for January 2026:**

KB5073723

**List of all patches/KB's validated on the AFFIRM PRONE System:**

KB5073723	KB5062557	KB5050008
KB5065225	KB5060531	KB5049608
KB5068791	KB5058392	KB5048661
KB5066143	KB5055519	KB5046615
KB5066586	KB5050533	KB5044277
KB5065428	KB890830	KB5044016
KB5063877	KB5053596	KB5043050
	KB5052000	KB5042215

KB5041578	KB5020627	KB5001342
KB5040430	KB5019966	KB5000859
KB5039879	KB5018419	KB5000822
KB5039217	KB5017315	KB5000685
KB5037765	KB5016884	KB4601393
KB5036896	KB5016623	KB4601345
KB5036604	KB5015811	KB4601060
KB5035849	KB5014692	KB4598480
KB5034768	KB5014553	KB4598230
KB5034127	KB5013941	KB4592440
KB5033904	KB5013641	KB4589208
KB5033371	KB5012647	KB4587735
KB5032196	KB5012170	KB4586793
KB5031984	KB5012128	KB4580325
KB5031361	KB5011503	KB4577668
KB5030214	KB5010351	KB4577667
KB5029931	KB5009557	KB4577586
KB5029376	KB5008873	KB4570333
KB5029247	KB5008218	KB4570332
KB5028168	KB5007206	KB4566424
KB5027222	KB5006672	KB4565349
KB5027131	KB5005568	KB4562562
KB5026362	KB5005112	KB4561608
KB5025229	KB5005030	KB4561600
KB5023702	KB5004244	KB4558998
KB5022840	KB5003711	KB4558997
KB5022511	KB5003646	KB4557397
KB5022286	KB5003243	KB4552924
KB5021237	KB5003171	KB4551853
KB5021126	KB5001407	KB4549949
KB5020866	KB5001404	KB4549947

KB4539571

KB4521862

KB4505225

KB4538461

KB4519338

KB4503308

KB4537759

KB4516115

KB4494174

KB4535680

KB4514366

KB2565063

KB4534273

KB4512937

KB2467173

KB4532691

KB4512578

KB4530715

KB4512577

KB4523205

KB4511553

KB4523204

KB4506998