



Date: February 10, 2026

Author: Software Engineering

System: Apex

Subsystem: Win Pro Embedded 10 LTSC 2019

Subject: Microsoft Win Pro Embedded 10 LTSC 2019 Patch release January 2026

Purpose:

To release the list of validated Microsoft critical patches on the Apex Win Pro Embedded 10 LTSC 2019 systems.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.21 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2019:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for January 2026:

KB5073723

List of all patches/KB's validated on the Apex System:

KB5073723	KB5066586	KB5055519
KB5071544	KB5065428	KB5054833
KB5065222	KB5063877	KB890830
KB5068791	KB5062557	KB5053596
KB5066143	KB5060531	KB5052000
	KB5058392	KB5050008

KB5049608	KB5027222	KB5008873
KB5048661	KB5027131	KB5008218
KB5046615	KB5026362	KB5007206
KB5044277	KB5025229	KB5006672
KB5044016	KB5023702	KB5005698
KB5043050	KB5022840	KB5005568
KB5042749	KB5022511	KB5005112
KB5041578	KB5022286	KB5005030
KB5040430	KB5021237	KB5003711
KB5039879	KB5021045	KB5002542
KB5039217	KB5020866	KB5002520
KB5037765	KB5019964	KB5002497
KB5036896	KB5018419	KB5002464
KB5036604	KB5017396	KB5002406
KB5035849	KB5017315	KB5002369
KB5034768	KB5016623	KB5002323
KB5034127	KB5015811	KB5002223
KB5033904	KB5014692	KB4589210
KB5033688	KB5014164	KB4589208
KB5033371	KB5013941	KB4583465
KB5032196	KB5013641	KB4583462
KB5031984	KB5012647	KB4577586
KB5031361	KB5012170	KB4535680
KB5030214	KB5012128	KB4494174
KB5029931	KB5011503	KB4132216
KB5029247	KB5010351	KB4022619
KB5028168	KB5009557	KB4018073