



**Date:** February 12, 2026

**Author:** Software Engineering

**System:** SELENIA DIMENSIONS / 3DIMENSIONS

**Subsystem:** Windows 10 LTSC 2019 (v1809) & Software Security Patches

**Subject:** Windows 10 LTSC 2019 (v1809) & Software Security Patch release  
January 2026

---

**Purpose:**

To release the list of validated Microsoft critical patches on the SELENIA DIMENSIONS / 3DIMENSIONS System.

**Latest .NET Versions Validated for the Product:**

**Note:** Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.21 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

- **.NET Core 2.1 / 2.2**

Hosting Bundle 2.1.30

**Latest SQLSERVER Versions Validated:**

**Note:** Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

**List of patches/KB's validated for January 2026:**

KB5066143

KB5073723

**List of all patches/KB's validated on the SELENIA DIMENSIONS / 3DIMENSIONS System:**

KB5066143	KB5073723	KB5071544
KB5065225	KB5043050	KB5028168
KB5068791	KB5042215	KB5027222
KB5066137	KB5041578	KB5027124
KB5066586	KB5040430	KB5026362
KB5065428	KB5039879	KB5025229
KB5063877	KB5039217	KB5023702
KB5062557	KB5037765	KB5022840
KB5060531	KB5036896	KB5022504
KB5058392	KB5036604	KB5022286
KB5055519	KB5035849	KB5021237
KB5050533	KB5034768	KB5021126
KB5049615	KB5034127	KB5020874
KB890830	KB5033904	KB5020615
KB5053596	KB5033371	KB5019966
KB5052000	KB5032196	KB5018419
KB5050008	KB5031990	KB5017315
KB5049608	KB5031361	KB5016884
KB5048661	KB5030214	KB5016623
KB5046615	KB5029925	KB5015811
KB5044277	KB5029376	KB5014692
KB5044016	KB5029247	KB5014553

KB5014164	KB5000685	KB4549947
KB5013941	KB4601393	KB4539571
KB5013641	KB4601345	KB4537759
KB5013626	KB4601060	KB4535680
KB5012647	KB4598480	KB4514366
KB5012170	KB4598230	KB4505225
KB5012119	KB4592440	KB4494174
KB5011503	KB4589208	KB4022619
KB5010351	KB4587735	KB4018073
KB5009557	KB4586793	KB2565063
KB5008873	KB4580325	KB2467173
KB5008218	KB4577668	
KB5007206	KB4577667	
KB5006672	KB4577586	
KB5005568	KB4570333	
KB5005112	KB4570332	
KB5005030	KB4566424	
KB5004244	KB4565349	
KB5003711	KB4562562	
KB5003646	KB4561608	
KB5003243	KB4561600	
KB5003171	KB4558998	
KB5001407	KB4558997	
KB5001404	KB4557397	
KB5001342	KB4552924	
KB5000859	KB4551853	
KB5000822	KB4549949	