



Date: February 12, 2026

Author: Software Engineering

System: SecurView Manager – Windows Server 2016

Subsystem: Windows 10 Server 2016 & Software Security Patches

Subject: Windows 10 Server 2016 & Software Security Patch release January 2026

Purpose:

To release the list of validated Microsoft critical patches on the SecurView Manager – Windows Server 2016 System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.21 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for January 2026:

KB5073447

KB5073722

List of all patches/KB's validated on the SecurView Manager – Windows Server 2016 System:

KB5073722	KB5073447	KB5071543
KB5068864	KB5066584	KB5065427
KB5070247	KB5066836	KB5065687
KB5063871	KB5062799	KB5060954
KB5062560	KB5061010	KB5058524

KB5055661	KB5040434	KB5030213
KB5055521	KB5039334	KB5029242
KB890830	KB5039214	KB5028169
KB5054006	KB5037763	KB5027219
KB5053594	KB5037016	KB5026363
KB5052006	KB5036899	KB5025228
KB5050109	KB5035962	KB5023788
KB5049993	KB5035855	KB5023697
KB5048671	KB5034862	KB5022838
KB5046612	KB5034767	KB5022289
KB5044293	KB5034119	KB5017396
KB5043124	KB5033373	KB5012170
KB5043051	KB5032391	KB5005698
KB5041773	KB5032197	KB4589210
KB5041576	KB5031362	KB4494175
KB5040562	KB5030504	KB4132216