



Date: February 12, 2026

Author: Software Engineering

System: Vision DXA

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release
January 2026

Purpose:

To release the list of validated Microsoft critical patches on the Vision DXA System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.21 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for January 2026:

KB5073723

List of all patches/KB's validated on the Vision DXA System:

KB5073723	KB5062557	KB5050008
KB5071544	KB5060531	KB5049608
KB5068791	KB5058392	KB5048661
KB5066143	KB5055519	KB5046615
KB5066586	KB890830	KB5044016
KB5065428	KB5053596	KB5041578
KB5063877	KB5052000	KB5040430

KB5039879	KB5030214	KB5020866
KB5039217	KB5029931	KB5020627
KB5037765	KB5029247	KB5019966
KB5036896	KB5028168	KB5018419
KB5036604	KB5027222	KB5017315
KB5035849	KB5027131	KB5016623
KB5034768	KB5026362	KB5015811
KB5034127	KB5025229	KB5013641
KB5033904	KB5023702	KB5012170
KB5033371	KB5022840	KB5005112
KB5032196	KB5022511	KB4577586
KB5031984	KB5022286	KB4577586
KB5031361	KB5021237	KB4535680