



Date: 08/06/2021

Author: Software Engineering

System: Apex

Subsystem: Microsoft Windows OS 10

Subject: Microsoft Windows 10 OS Patch release July 2021

Purpose:

To release the list of validated Microsoft critical patches on the Apex systems.

Note 1: If you are using the following versions of APEX, if you are required to update to the latest Hologic validated Windows updates, you must ensure your system(s) have been patched with the “Sentinel Run-Time 7.80.0.0” which can be located in the Windows Control Panel “Program and Features”. If needed you can install the patch by ordering “APEX Service Dongle Patch Release” (ASY-12610). Please contact Hologic Support if you have any questions.

System Type	Software Versions
Horizon	APEX v5.6.0.7
Workstation	APEX v5.6.0.7

Note 2: Below listed patches for this month are cumulative and include “Print Nightmare” fixes. Additional testing to test “Print Nightmare” fixes has been performed for this particular patch validation.

List of patches/KB’s validated for July 2021:

KB5004238

List of all patches/KB’s validated on the Apex System:

KB890830	KB4601392	KB4580325
KB5004238	KB4601318	KB4577586
KB5003638	KB4598243	KB4577015
KB5003197	KB4593226	KB4576750
KB5001402	KB4589210	KB4571694
KB5001347	KB4586830	KB4565912
KB5001078	KB4583462	KB4565511
KB5000803	KB4580346	KB4562561

KB4561616

KB4561600

KB4556813

KB4550994

KB4550929

KB4540723

KB4540670

KB4537764

KB4537759

KB4535680

KB4535288

KB4534271

KB4530689

KB4520724

KB4516115

KB4516044

KB4512574

KB4505422

KB4503537

KB4503308

KB4503267

KB4498947

KB4494175

KB4346087

KB4132216

KB4022619

KB4018073

KB2565063

KB2467173