



**Date:** August 06, 2021

**Author:** Software Engineering

**System:** SELENIA DIMENSIONS / 3DIMENSIONS

**Subsystem:** Windows 10 LTSC 2019 (v1809) & Software Security Patches

**Subject:** Windows 10 LTSC 2019 (v1809) & Software Security Patch release July 2021

**Purpose:**

To release the list of validated Microsoft critical patches on the SELENIA DIMENSIONS / 3DIMENSIONS System.

Note: Below listed patches for this month are cumulative and include “Print Nightmare” fixes. Additional testing to test “Print Nightmare” fixes has been performed for this particular patch validation.

**List of patches/KB’s validated for July 2021:**

KB5004244

**List of all patches/KB’s validated on the SELENIA DIMENSIONS / 3DIMENSIONS System:**

KB890830	KB5000859	KB4589208
KB5004244	KB5000822	KB4587735
KB5003711	KB5000685	KB4586793
KB5003646	KB4601393	KB4580325
KB5003243	KB4601345	KB4577668
KB5003171	KB4601060	KB4577667
KB5001407	KB4598480	KB4577586
KB5001404	KB4598230	KB4570333
KB5001342	KB4592440	KB4570332

KB4566424

KB4557397

KB4535680

KB4565349

KB4552924

KB4514366

KB4562562

KB4551853

KB4505225

KB4561608

KB4549949

KB4494174

KB4561600

KB4549947

KB2565063

KB4558998

KB4539571

KB2467173

KB4558997

KB4537759