



Date: August 06, 2021

Author: Software Engineering

System: TRIDENT

Subsystem: Windows 7 & Software Security Patches

Subject: Windows 7 & Software Security Patch release July 2021

Purpose:

To release the list of validated Microsoft critical patches on the TRIDENT System.

Note 1:

- Please check if January 2020 Windows 7 Monthly Patch Disk is installed on system, if not then first Install January 2020 Windows 7 Monthly Patch Disk (SFW-03075).
- Apply SFW-03338 Rev. 001 (Windows 7 ESU License Installation) to the system prior to installation of November 2020 Windows 7 Monthly Patch Disk (SFW-03339) and April 2021 Windows 7 Monthly Patch Disk (SFW-03430) to allow for updates to be installed, due to end of support for Microsoft Windows 7.

Note 2: Below listed patches for this month are cumulative and include “Print Nightmare” fixes. Additional testing to test “Print Nightmare” fixes has been performed for this particular patch validation.

List of patches/KB’s validated for July 2021:

| | | |
|-----------|-----------|-----------|
| KB5004951 | KB5004307 | KB5004120 |
| KB5004378 | KB5004289 | |

List of all patches/KB’s validated on the TRIDENT System:

| | | |
|----------|-----------|-----------|
| KB976932 | KB5004951 | KB5004378 |
|----------|-----------|-----------|

| | | |
|-----------|-----------|-----------|
| KB5004307 | KB4578963 | KB4538483 |
| KB5004289 | KB4578952 | KB4537829 |
| KB5004120 | KB4577053 | KB4537813 |
| KB5003694 | KB4577051 | KB4536952 |
| KB5003667 | KB4576612 | KB4534314 |
| KB5003233 | KB4575903 | KB4534310 |
| KB5003228 | KB4571729 | KB4533012 |
| KB5001848 | KB4571719 | KB4532945 |
| KB5001392 | KB4570673 | KB4532932 |
| KB5001335 | KB4569775 | KB4531786 |
| KB5000851 | KB4569767 | KB4530734 |
| KB5000841 | KB4565623 | KB4530692 |
| KB4601363 | KB4565612 | KB4525235 |
| KB4601347 | KB4565539 | KB4525233 |
| KB4600945 | KB4565524 | KB4523206 |
| KB4598289 | KB4565354 | KB4523206 |
| KB4598279 | KB4562030 | KB4520003 |
| KB4597239 | KB4561669 | KB4519976 |
| KB4592510 | KB4561643 | KB4516655 |
| KB4592503 | KB4556843 | KB4516065 |
| KB4592471 | KB4556836 | KB4516033 |
| KB4586827 | KB4556798 | KB4515854 |
| KB4586805 | KB4555449 | KB4512506 |
| KB4585204 | KB4552940 | KB4512486 |
| KB4580970 | KB4552919 | KB4511516 |
| KB4580387 | KB4550965 | KB4507456 |
| KB4580345 | KB4541500 | KB4507449 |

| | | |
|-----------|-----------|-----------|
| KB4507004 | KB4462915 | KB4012212 |
| KB4506997 | KB4457145 | KB3212642 |
| KB4505422 | KB4457044 | KB3205394 |
| KB4503292 | KB4457035 | KB3197867 |
| KB4503269 | KB4343899 | KB3192391 |
| KB4499175 | KB4339284 | KB3177467 |
| KB4499164 | KB4338823 | KB3172605 |
| KB4495606 | KB4284867 | KB3161949 |
| KB4495588 | KB4103712 | KB3159398 |
| KB4493472 | KB4096418 | KB3156016 |
| KB4493448 | KB4093108 | KB3150220 |
| KB4490628 | KB4074837 | KB3139914 |
| KB4489885 | KB4074587 | KB3139398 |
| KB4489878 | KB4057114 | KB3138910 |
| KB4486564 | KB4056897 | KB3138612 |
| KB4486563 | KB4054521 | KB3133977 |
| KB4483458 | KB4048960 | KB3127220 |
| KB4483451 | KB4041678 | KB3126587 |
| KB4483187 | KB4038779 | KB3124275 |
| KB4480970 | KB4034679 | KB3122648 |
| KB4480960 | KB4025337 | KB3110329 |
| KB4480063 | KB4022722 | KB3109560 |
| KB4480055 | KB4022619 | KB3109103 |
| KB4474419 | KB4020322 | KB3108664 |
| KB4471328 | KB4019990 | KB3108371 |
| KB4468323 | KB4019263 | KB3101722 |
| KB4467106 | KB4015546 | KB3097989 |

| | | |
|-----------|-----------|-----------|
| KB3092601 | KB2992611 | KB2862330 |
| KB3086255 | KB2991963 | KB2862152 |
| KB3080079 | KB2984972 | KB2861698 |
| KB3078601 | KB2979596 | KB2840631 |
| KB3075220 | KB2978742 | KB2834140 |
| KB3074543 | KB2977292 | KB2813430 |
| KB3072305 | KB2973201 | KB2813347 |
| KB3071756 | KB2973112 | KB2807986 |
| KB3067903 | KB2972211 | KB2789645 |
| KB3061518 | KB2972100 | KB2786081 |
| KB3060716 | KB2968294 | KB2770660 |
| KB3059317 | KB2943357 | KB2758857 |
| KB3055642 | KB2937610 | KB2742599 |
| KB3046269 | KB2931356 | KB2736422 |
| KB3046017 | KB2912390 | KB2731771 |
| KB3045685 | KB2911501 | KB2729452 |
| KB3042058 | KB2900986 | KB2729094 |
| KB3037574 | KB2894854 | KB2727528 |
| KB3035126 | KB2894844 | KB2719033 |
| KB3030377 | KB2893294 | KB2716513 |
| KB3023215 | KB2888049 | KB2705219 |
| KB3022777 | KB2884256 | KB2698365 |
| KB3021674 | KB2882822 | KB2690533 |
| KB3019978 | KB2871997 | KB2685939 |
| KB3018238 | KB2868038 | KB2670838 |
| KB3011780 | KB2864202 | KB2667402 |
| KB3010788 | KB2862335 | KB2656356 |

KB2654428

KB2653956

KB2639308

KB2631813

KB2621440

KB2620704

KB2604115

KB2585542

KB2579686

KB2565063

KB2564958

KB2560656

KB2552343

KB2533623

KB2533552

KB2532531

KB2510531

KB2506212

KB2506014

KB2491683

KB2479943

KB2467173