



Date: August 8, 2022

Author: Software Engineering

System: SecurView

Subsystem: Windows 10 LTSC 2016 (v1607) & Software Security Patches

Subject: Windows 10 LTSC 2016 (v1607) & Software Security Patch release July 2022

Purpose:

To release the list of validated Microsoft critical patches on the SecurView System.

List of patches/KB's validated for July 2022:

KB5016058 KB5015808 KB4577586

List of all patches/KB's validated on the SecurView System:

KB890830	KB5008877	KB5000803
KB5016058	KB5008207	KB5000685
KB5015808	KB5007192	KB4601392
KB5014702	KB5006669	KB4601318
KB5014026	KB5005698	KB4601051
KB5013952	KB5005573	KB4598243
KB5013625	KB5005043	KB4593226
KB5012596	KB5004238	KB4589210
KB5012118	KB5003638	KB4586830
KB5011570	KB5003197	KB4580346
KB5011495	KB5001402	KB4578969
KB5010359	KB5001347	KB4577586
KB5009546	KB5001078	KB4577015

KB4576750	KB4503267
KB4571694	KB4498947
KB4565912	KB4494440
KB4565511	KB4494175
KB4562561	KB4493472
KB4561616	KB4493470
KB4557397	KB4493448
KB4556813	KB4493435
KB4550994	KB4489882
KB4550929	KB4485447
KB4540723	KB4346087
KB4540670	KB4132216
KB4537764	KB4091664
KB4535680	KB2565063
KB4534271	KB2467173
KB4530689	
KB4525236	
KB4521858	
KB4520724	
KB4519998	
KB4516044	
KB4512574	
KB4512517	
KB4509091	
KB4507460	
KB4505225	
KB4503537	