



Date: August 07, 2023

Author: Software Engineering

System: Apex

Subsystem: Microsoft Win Pro Embedded 10 LTSP 2016

Subject: Microsoft Win Pro Embedded 10 LTSP 2016 Patch release July 2023

Purpose:

To release the list of validated Microsoft critical patches on the Apex Win Pro Embedded 10 LTSP 2016 systems.

Note 1: If you are using the following versions of APEX, if you are required to update to the latest Hologic validated Windows updates, you must ensure your system(s) have been patched with the “Sentinel Run-Time 7.80.0.0” which can be located in the Windows Control Panel “Program and Features”. If needed you can install the patch by ordering “APEX Service Dongle Patch Release” (ASY-12610). Please contact Hologic Support if you have any questions.

System Type	Software Versions
Horizon	APEX v5.6.0.7
Workstation	APEX v5.6.0.7

List of patches/KB’s validated for July 2023:

KB5028169 KB5002406

List of all patches/KB’s validated on the Apex System:

- | | | |
|-----------|-----------|-----------|
| KB890830 | KB5021045 | KB5016058 |
| KB5028169 | KB5020627 | KB5015808 |
| KB5027219 | KB5019966 | KB5014702 |
| KB5026363 | KB5019964 | KB5014164 |
| KB5025228 | KB5018411 | KB5014026 |
| KB5023788 | KB5017396 | KB5013952 |
| KB5023697 | KB5017305 | KB5012596 |
| KB5022838 | KB5017095 | KB5012170 |
| KB5022289 | KB5016884 | KB5011570 |
| KB5021235 | KB5016622 | KB5011495 |

KB5010359	KB4601318	KB4537764
KB5009557	KB4598243	KB4537759
KB5008873	KB4593226	KB4535680
KB5008207	KB4589210	KB4535288
KB5007192	KB4589208	KB4534271
KB5006669	KB4586830	KB4530689
KB5005698	KB4583465	KB4520724
KB5005573	KB4583462	KB4516115
KB5005112	KB4580346	KB4516044
KB5005043	KB4580325	KB4512574
KB5004238	KB4577586	KB4505422
KB5003638	KB4577015	KB4503537
KB5003197	KB4576750	KB4503308
KB5002406	KB4571694	KB4503267
KB5002369	KB4565912	KB4498947
KB5002323	KB4565511	KB4494175
KB5002223	KB4562561	KB4346087
KB5001402	KB4561616	KB4132216
KB5001347	KB4561600	KB4052908
KB5001078	KB4556813	KB4022619
KB5000803	KB4550994	KB4018073
KB5000685	KB4550929	KB2565063
KB5000645	KB4540723	KB2467173
KB4601392	KB4540670	