



Date: August 11, 2025

Author: Software Engineering

System: Vision DXA

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release July 2025

Purpose:

To release the list of validated Microsoft critical patches on the Vision DXA System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.14 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for July 2025:

KB5062557

List of all patches/KB's validated on the Vision DXA System:

KB5062557	KB5049608	KB5037765
KB5060531	KB5048661	KB5036896
KB5058392	KB5046615	KB5036604
KB5055519	KB5044016	KB5035849
KB890830	KB5041578	KB5034768
KB5053596	KB5040430	KB5034127
KB5052000	KB5039879	KB5033904
KB5050008	KB5039217	KB5033371

KB5032196	KB5025229	KB5017315
KB5031984	KB5023702	KB5016623
KB5031361	KB5022840	KB5015811
KB5030214	KB5022511	KB5013641
KB5029931	KB5022286	KB5012170
KB5029247	KB5021237	KB5005112
KB5028168	KB5020866	KB4589208
KB5027222	KB5020627	KB4577586
KB5027131	KB5019966	KB4535680
KB5026362	KB5018419	