



The Science of Sure

Date: *June 22, 2020*

Author: Software Engineering

System: SELENIA DIMENSIONS / 3DIMENSIONS

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release June 2020

Purpose:

To release the list of validated Microsoft critical patches on the SELENIA DIMENSIONS / 3DIMENSIONS System.

List of patches/KB's validated for June 2020:

| | | |
|-----------|-----------|-----------|
| KB4562562 | KB4561608 | KB4561600 |
|-----------|-----------|-----------|

List of all patches/KB's validated on the SELENIA DIMENSIONS / 3DIMENSIONS System:

| | | |
|-----------|-----------|-----------|
| KB4562562 | KB4539571 | KB4494174 |
| KB4561608 | KB4514366 | KB2565063 |
| KB4561600 | KB4505225 | KB2467173 |