



Date: July 18, 2023

Author: Software Engineering

System: InSight

Subsystem: Microsoft Windows 10 Enterprise 2019 LTSC

Subject: Microsoft Windows 10 Enterprise 2019 LTSC Patch release June 2023

Purpose:

To release the list of validated Microsoft critical patches on the InSight systems.

System Type	Software Versions
Workstation	InSight v6.1.4

List of patches/KB's validated for June 2023:

KB5027222 KB5027131

List of all patches/KB's validated on the Insight System:

KB890830	KB5023702	KB5004424
KB5027222	KB5022840	KB4589208
KB5027131	KB5022511	KB4577586
KB5026362	KB5012170	
KB5025229	KB5005112	