

Date: July 16, 2025

Author: Software Engineering

System: AFFIRM PRONE

Subsystem: Windows 10 & Software Security Patches

Subject: Windows 10 & Software Security Patch release June 2025

Purpose:

To release the list of validated Microsoft critical patches on the AFFIRM PRONE System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.14 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for June 2025:

KB5060531

List of all patches/KB's validated on the AFFIRM PRONE System:

KB5060531	KB5044016	KB5034127
KB5058392	KB5043050	KB5033904
KB5055519	KB5042215	KB5033371
KB5050533	KB5041578	KB5032196
KB890830	KB5040430	KB5031984
KB5053596	KB5039879	KB5031361
KB5052000	KB5039217	KB5030214
KB5050008	KB5037765	KB5029931
KB5049608	KB5036896	KB5029376
KB5048661	KB5036604	KB5029247
KB5046615	KB5035849	KB5028168
KB5044277	KB5034768	KB5027222

KB5027131	KB5005568	KB4562562
KB5026362	KB5005112	KB4561608
KB5025229	KB5005030	KB4561600
KB5023702	KB5004244	KB4558998
KB5022840	KB5003711	KB4558997
KB5022511	KB5003646	KB4557397
KB5022286	KB5003243	KB4552924
KB5021237	KB5003171	KB4551853
KB5021126	KB5001407	KB4549949
KB5020866	KB5001404	KB4549947
KB5020627	KB5001342	KB4539571
KB5019966	KB5000859	KB4538461
KB5018419	KB5000822	KB4537759
KB5017315	KB5000685	KB4535680
KB5016884	KB4601393	KB4534273
KB5016623	KB4601345	KB4532691
KB5015811	KB4601060	KB4530715
KB5014692	KB4598480	KB4523205
KB5014553	KB4598230	KB4523204
KB5013941	KB4592440	KB4521862
KB5013641	KB4589208	KB4519338
KB5012647	KB4587735	KB4516115
KB5012170	KB4586793	KB4514366
KB5012128	KB4580325	KB4512937
KB5011503	KB4577668	KB4512578
KB5010351	KB4577667	KB4512577
KB5009557	KB4577586	KB4511553
KB5008873	KB4570333	KB4506998
KB5008218	KB4570332	KB4505225
KB5007206	KB4566424	KB4503308
KB5006672	KB4565349	KB4494174

KB2565063

KB2467173