

Date: July 16, 2025

Author: Software Engineering

System: AWM

Subsystem: Windows Server 2016 & Software Security Patches

Subject: Windows Server 2016 & Software Security Patch release June 2025

Purpose:

To release the list of validated Microsoft critical patches on the AWM System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.14 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

- **.NET Core 2.1 / 2.2**

Hosting Bundle 2.1.30

List of patches/KB's validated for June 2025:

KB5061010 KB5060954

List of all patches/KB's validated on the AWM System:

KB5060954	KB5050109	KB5040562
KB5061010	KB5049993	KB5040434
KB5058392	KB5048671	KB5039879
KB5049608	KB5046612	KB5039334
KB5055661	KB5044293	KB5039214
KB5055521	KB5043124	KB5037763
KB5050533	KB5043051	KB5037016
KB890830	KB5042215	KB5036899
KB5054006	KB5041773	KB5035962
KB5053594	KB5041578	KB5035855
KB5052006	KB5041576	KB5034862

KB5034767	KB5014553	KB4589210
KB5034119	KB5014026	KB4589208
KB5033373	KB5013952	KB4586830
KB5032391	KB5012596	KB4580346
KB5032197	KB5012170	KB4577586
KB5031362	KB5011570	KB4577015
KB5030504	KB5011495	KB4576750
KB5030213	KB5010359	KB4571694
KB5029376	KB5009557	KB4565912
KB5029242	KB5009546	KB4565511
KB5028169	KB5008873	KB4562561
KB5027219	KB5008207	KB4561616
KB5026363	KB5007192	KB4557397
KB5025228	KB5006669	KB4551853
KB5023788	KB5005698	KB4550929
KB5023697	KB5005573	KB4549947
KB5022838	KB5005112	KB4540723
KB5022289	KB5005043	KB4539571
KB5021235	KB5004238	KB4537759
KB5021126	KB5003638	KB4535680
KB5019964	KB5003197	KB4514366
KB5018411	KB5001402	KB4505225
KB5017396	KB5001347	KB4494175
KB5017305	KB5001078	KB4494174
KB5017095	KB5000803	KB4132216
KB5016884	KB5000685	KB4132214
KB5016622	KB4601392	KB2565063
KB5016058	KB4601318	KB2467173
KB5015808	KB4598243	
KB5014702	KB4593226	