

Date: July 16, 2025

Author: Software Engineering

System: Apex

Subsystem: Microsoft Win Pro Embedded 10 LTSC 2016

Subject: Microsoft Win Pro Embedded 10 LTSC 2016 Patch release June 2025

Purpose:

To release the list of validated Microsoft critical patches on the Apex Win Pro Embedded 10 LTSC 2016 systems.

Note: If you are using the following versions of APEX, if you are required to update to the latest Hologic validated Windows updates, you must ensure your system(s) have been patched with the "Sentinel Run-Time 7.80.0.0" which can be located in the Windows Control Panel "Program and Features". If needed you can install the patch by ordering "APEX Service Dongle Patch Release" (ASY-12610). Please contact Hologic Support if you have any questions.

System Type	Software Versions
Horizon	APEX v5.6.0.7
Workstation	APEX v5.6.0.7

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.14 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for June 2025:

KB5061010

KB5060954

List of all patches/KB's validated on the Apex System:

KB5061010

KB5055521

KB5052006

KB5060954

KB5054833

KB5050109

KB5058524

KB890830

KB5049993

KB5058383

KB5054006

KB5048671

KB5055661

KB5053594

KB5046612

KB5044293	KB5030213	KB5018411
KB5043124	KB5029931	KB5017396
KB5043051	KB5029247	KB5017315
KB5042749	KB5029242	KB5017305
KB5041773	KB5029185	KB5017095
KB5041576	KB5028169	KB5016884
KB5040562	KB5028168	KB5016623
KB5040434	KB5027222	KB5016622
KB5039334	KB5027219	KB5016058
KB5039214	KB5027131	KB5015811
KB5037763	KB5026363	KB5015808
KB5037016	KB5026362	KB5014702
KB5036899	KB5025229	KB5014692
KB5035962	KB5025228	KB5014553
KB5035855	KB5023788	KB5014026
KB5034862	KB5023702	KB5013952
KB5034767	KB5023697	KB5013941
KB5034127	KB5022840	KB5013641
KB5034119	KB5022838	KB5012647
KB5033904	KB5022511	KB5012596
KB5033688	KB5022289	KB5012170
KB5033373	KB5022286	KB5012128
KB5033371	KB5021237	KB5011570
KB5032391	KB5021235	KB5011503
KB5032197	KB5021126	KB5011495
KB5032196	KB5021045	KB5010359
KB5031984	KB5020866	KB5010351
KB5031362	KB5020627	KB5009557
KB5031361	KB5019966	KB5009546
KB5030504	KB5019964	KB5008873
KB5030214	KB5018419	KB5008218

KB5008207	KB5000822	KB4565511
KB5007206	KB5000803	KB4565349
KB5007192	KB5000685	KB4562562
KB5006672	KB4601393	KB4562561
KB5006669	KB4601392	KB4561616
KB5005698	KB4601345	KB4561608
KB5005573	KB4601318	KB4561600
KB5005568	KB4601060	KB4558998
KB5005112	KB4598480	KB4558997
KB5005043	KB4598243	KB4557397
KB5005030	KB4598230	KB4556813
KB5004244	KB4593226	KB4552924
KB5004238	KB4592440	KB4551853
KB5003711	KB4589210	KB4550994
KB5003646	KB4589208	KB4550929
KB5003638	KB4587735	KB4549949
KB5003243	KB4586830	KB4549947
KB5003197	KB4586793	KB4540723
KB5003171	KB4583465	KB4540670
KB5002542	KB4580346	KB4539571
KB5002520	KB4580325	KB4538461
KB5002497	KB4577668	KB4537764
KB5002464	KB4577667	KB4537759
KB5002406	KB4577586	KB4535680
KB5001407	KB4577015	KB4534273
KB5001404	KB4576750	KB4534271
KB5001402	KB4571694	KB4532691
KB5001347	KB4570333	KB4530715
KB5001342	KB4570332	KB4530689
KB5001078	KB4566424	KB4525236
KB5000859	KB4565912	KB4523205

KB4523204	KB4512574	KB4493470
KB4521862	KB4512517	KB4489882
KB4521858	KB4511553	KB4487026
KB4520724	KB4509091	KB4485447
KB4519998	KB4507460	KB4346087
KB4519338	KB4506998	KB4132216
KB4516115	KB4505225	KB4091664
KB4516044	KB4503308	KB4022619
KB4514366	KB4498947	KB4018073
KB4512937	KB4494440	KB2565063
KB4512578	KB4494175	KB2467173
KB4512577	KB4494174	