



Date: July 16, 2025

Author: Software Engineering

System: Cenova

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release June 2025

Purpose:

To release the list of validated Microsoft critical patches on the Cenova System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.14 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for June 2025:

KB5060531

List of all patches/KB's validated on the Cenova System:

KB5060531	KB5049608	KB5040430
KB5058392	KB5048661	KB5039879
KB5055519	KB5046615	KB5039217
KB5050533	KB5044277	KB5037765
KB890830	KB5044016	KB5036896
KB5053596	KB5043050	KB5036604
KB5052000	KB5042215	KB5035849
KB5050008	KB5041578	KB5034768

KB5034127	KB5019966	KB5003711
KB5033904	KB5018419	KB5003646
KB5033371	KB5017315	KB5003638
KB5032196	KB5016884	KB5003243
KB5031984	KB5016623	KB5003171
KB5031361	KB5015811	KB5001407
KB5030214	KB5014692	KB5001404
KB5029931	KB5014553	KB5001402
KB5029376	KB5013941	KB5001342
KB5029247	KB5013641	KB5000859
KB5028168	KB5012647	KB5000822
KB5027222	KB5012170	KB5000685
KB5027131	KB5012128	KB4601393
KB5027124	KB5011503	KB4601345
KB5026362	KB5010351	KB4601060
KB5025229	KB5009557	KB4589210
KB5023702	KB5008873	KB4589208
KB5022840	KB5008218	KB4580325
KB5022511	KB5007206	KB4577586
KB5022286	KB5006672	KB4557397
KB5021237	KB5005568	KB4539571
KB5021126	KB5005112	KB4535680
KB5020866	KB5005030	KB4494174
KB5020627	KB5004244	KB4465065