



**Date:** July 09, 2025

**Author:** Software Engineering

**System:** Path Vision

**Subsystem:** Windows 10 LTSC 2019 (v1809) & Software Security Patches

**Subject:** Windows 10 LTSC 2019 (v1809) & Software Security Patch release June 2025

---

---

**Purpose:**

To release the list of validated Microsoft critical patches on the Path Vision System.

Note: Core Vision and BioVision are on similar product line as Path Vision and both use the same executable so patches validated on Path Vision system can be installed on Core Visions system and BioVision as well.

**Latest .NET Versions Validated for the Product:**

**Note:** Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.14 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

**List of patches/KB's validated for June 2025:**

KB5060531

**List of all patches/KB's validated on the Path Vision System:**

KB5060531	KB5052000	KB5044277
KB5058392	KB5050008	KB5044016
KB5055519	KB5049608	KB5043050
KB890830	KB5048661	KB5041578
KB5053596	KB5046615	KB5040430

KB5039879	KB5029931	KB5018419
KB5039217	KB5029247	KB5017315
KB5037765	KB5028168	KB5016623
KB5036896	KB5027222	KB5015811
KB5036604	KB5027131	KB5014692
KB5035849	KB5026362	KB5013941
KB5034768	KB5025229	KB5013641
KB5034127	KB5023702	KB5012170
KB5033904	KB5022840	KB5005112
KB5033371	KB5022511	KB4589208
KB5032196	KB5021237	KB4577586
KB5031984	KB5020866	KB4535680
KB5031361	KB5020627	KB4494174
KB5030214	KB5019966	