



Date: July 16, 2025

Author: Software Engineering

System: SELENIA DIMENSIONS / 3DIMENSIONS

Subsystem: Windows 10 LTSB 2016 (v1607) & Software Security Patches

Subject: Windows 10 LTSB 2016 (v1607) & Software Security Patch release June 2025

Purpose:

To release the list of validated Microsoft critical patches on the SELENIA DIMENSIONS / 3DIMENSIONS System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.14 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

- **.NET Core 2.1 / 2.2**

Hosting Bundle 2.1.30

List of patches/KB's validated for June 2025:

KB5061010

KB5060954

List of all patches/KB's validated on the SELENIA DIMENSIONS / 3DIMENSIONS System:

KB5060954

KB5058383

KB5050533

KB5061010

KB5055661

KB890830

KB5058524

KB5055521

KB5054006

KB5053594	KB5031362	KB5013952
KB5052006	KB5030504	KB5013625
KB5050109	KB5030213	KB5012596
KB5049993	KB5029376	KB5012170
KB5048671	KB5029242	KB5012118
KB5046612	KB5028169	KB5011570
KB5044293	KB5027219	KB5011495
KB5043124	KB5026363	KB5010359
KB5043051	KB5025228	KB5009546
KB5042215	KB5023788	KB5008877
KB5041773	KB5023697	KB5008207
KB5041576	KB5022838	KB5007192
KB5040562	KB5022289	KB5006669
KB5040434	KB5021235	KB5005698
KB5039334	KB5021126	KB5005573
KB5039214	KB5019964	KB5005043
KB5037763	KB5018411	KB5004238
KB5037016	KB5017396	KB5003638
KB5036899	KB5017305	KB5003197
KB5035962	KB5017095	KB5001402
KB5035855	KB5016884	KB5001347
KB5034862	KB5016622	KB5001078
KB5034767	KB5016058	KB5000803
KB5034119	KB5015808	KB5000685
KB5033373	KB5014702	KB4601392
KB5032391	KB5014553	KB4601318
KB5032197	KB5014026	KB4601051

KB4598243	KB4519998
KB4593226	KB4516044
KB4589210	KB4512574
KB4586830	KB4512517
KB4580346	KB4509091
KB4578969	KB4507460
KB4577586	KB4505225
KB4577015	KB4503537
KB4576750	KB4503267
KB4571694	KB4498947
KB4565912	KB4494440
KB4565511	KB4494175
KB4562561	KB4493472
KB4561616	KB4493470
KB4557397	KB4493448
KB4556813	KB4493435
KB4550994	KB4489882
KB4550929	KB4485447
KB4540723	KB4346087
KB4540670	KB4132216
KB4537764	KB4091664
KB4535680	KB2565063
KB4534271	KB2467173
KB4530689	
KB4525236	
KB4521858	
KB4520724	