



Date: July 16, 2025

Author: Software Engineering

System: SELENIA DIMENSIONS / 3DIMENSIONS

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release June 2025

Purpose:

To release the list of validated Microsoft critical patches on the SELENIA DIMENSIONS / 3DIMENSIONS System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.14 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

- **.NET Core 2.1 / 2.2**

Hosting Bundle 2.1.30

List of patches/KB's validated for June 2025:

KB5060531

List of all patches/KB's validated on the SELENIA DIMENSIONS / 3DIMENSIONS System:

KB5060531	KB5049615	KB5050008
KB5058392	KB890830	KB5049608
KB5055519	KB5053596	KB5048661
KB5050533	KB5052000	KB5046615

KB5044277	KB5025229	KB5008873
KB5044016	KB5023702	KB5008218
KB5043050	KB5022840	KB5007206
KB5042215	KB5022504	KB5006672
KB5041578	KB5022286	KB5005568
KB5040430	KB5021237	KB5005112
KB5039879	KB5021126	KB5005030
KB5039217	KB5020874	KB5004244
KB5037765	KB5020615	KB5003711
KB5036896	KB5019966	KB5003646
KB5036604	KB5018419	KB5003243
KB5035849	KB5017315	KB5003171
KB5034768	KB5016884	KB5001407
KB5034127	KB5016623	KB5001404
KB5033904	KB5015811	KB5001342
KB5033371	KB5014692	KB5000859
KB5032196	KB5014553	KB5000822
KB5031990	KB5014164	KB5000685
KB5031361	KB5013941	KB4601393
KB5030214	KB5013641	KB4601345
KB5029925	KB5013626	KB4601060
KB5029376	KB5012647	KB4598480
KB5029247	KB5012170	KB4598230
KB5028168	KB5012119	KB4592440
KB5027222	KB5011503	KB4589208
KB5027124	KB5010351	KB4587735
KB5026362	KB5009557	KB4586793

KB4580325

KB2467173

KB4577668

KB4577667

KB4577586

KB4570333

KB4570332

KB4566424

KB4565349

KB4562562

KB4561608

KB4561600

KB4558998

KB4558997

KB4557397

KB4552924

KB4551853

KB4549949

KB4549947

KB4539571

KB4537759

KB4535680

KB4514366

KB4505225

KB4494174

KB4022619

KB4018073

KB2565063