



Date: July 16, 2025

Author: Software Engineering

System: SecurView Manager – Windows Server 2016

Subsystem: Windows 10 Server 2016 & Software Security Patches

Subject: Windows 10 Server 2016 & Software Security Patch release June 2025

Purpose:

To release the list of validated Microsoft critical patches on the SecurView Manager – Windows Server 2016 System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.14 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for June 2025:

KB5061010 KB5060954

List of all patches/KB's validated on the SecurView Manager – Windows Server 2016 System:

KB5061010	KB5054006	KB5046612
KB5060954	KB5053594	KB5044293
KB5058524	KB5052006	KB5043124
KB5055661	KB5050109	KB5043051
KB5055521	KB5049993	KB5041773
KB890830	KB5048671	KB5041576

KB5040562	KB5034119	KB5025228
KB5040434	KB5033373	KB5023788
KB5039334	KB5032391	KB5023697
KB5039214	KB5032197	KB5022838
KB5037763	KB5031362	KB5022289
KB5037016	KB5030504	KB5017396
KB5036899	KB5030213	KB5012170
KB5035962	KB5029242	KB5005698
KB5035855	KB5028169	KB4589210
KB5034862	KB5027219	KB4494175
KB5034767	KB5026363	KB4132216