



Date: *March 30, 2020*

Author: Software Engineering

System: SELENIA DIMENSIONS / 3DIMENSIONS

Subsystem: Windows 10 & Software Security Patches

Subject: Windows 10 & Software Security Patch release March 2020

Purpose:

To release the list of validated Microsoft critical patches on the SELENIA DIMENSIONS / 3DIMENSIONS System.

List of patches/KB's validated for March 2020:

KB4540723 KB4540670

List of all patches/KB's validated on the SELENIA DIMENSIONS / 3DIMENSIONS System:

KB4540723	KB4505225	KB2565063
KB4540670	KB4494175	KB2467173
KB4520724	KB4132216	
KB2467173		