



Date: 03/22/2021

Author: Software Engineering

System: Apex

Subsystem: Microsoft Windows OS 7

Subject: Microsoft Windows 7 OS Patch release January 2021

Purpose:

To release the list of validated Microsoft critical patches on the Apex systems.

Note 1: If you are using the following versions of APEX, if you are required to update to the latest Hologic validated Windows updates, you must ensure your system(s) have been patched with the “Sentinel Run-Time 7.80.0.0” which can be located in the Windows Control Panel “Program and Features”. If needed you can install the patch by ordering “APEX Service Dongle Patch Release” (ASY-12610). Please contact Hologic Support if you have any questions.

Note 2:

- Please check if January 2020 Windows 7 Monthly Patch Disk is installed on system, if not then first Install January 2020 Windows 7 Monthly Patch Disk (SFW-03075).
- Apply SFW-03338 Rev. 001 (Windows 7 ESU License Installation) to the system prior to installation of January 2021 Windows 7 Monthly Patch Disk (SFW-03375) to allow for updates to be installed, due to end of support for Microsoft Windows 7.
- For now only the January 2021 Win7 monthly patch disk (SFW-03375) should be applied to Apex Windows 7 systems. February 2021 monthly patch disk (SFW-03390) should not be applied on APEX Windows 7 systems, until further instructions.

System Type	Software Versions
Horizon	APEX v5.6.0.7
Workstation	APEX v5.6.0.7

List of patches/KB’s validated for January 2021:

KB4598289

KB4598279

KB4597239

List of all patches/KB's validated on the Apex System:

KB976932	KB4534310	KB4457044
KB4598289	KB4534251	KB4457035
KB4598279	KB4532945	KB4343899
KB4597239	KB4532932	KB4338823
KB4592510	KB4530692	KB4284867
KB4592503	KB4525233	KB4103712
KB4592471	KB4520003	KB4096418
KB4586827	KB4519108	KB4093108
KB4586805	KB4516033	KB4074587
KB4586768	KB4512486	KB4057116
KB4585204	KB4507704	KB4056897
KB4580970	KB4507456	KB4054521
KB4580387	KB4507004	KB4048960
KB4578952	KB4503269	KB4041678
KB4578623	KB4501226	KB4038779
KB4577053	KB4499175	KB4034679
KB4575903	KB4493448	KB4025337
KB4571719	KB4490628	KB4022722
KB4566371	KB4490128	KB4022619
KB4565539	KB4489885	KB4019990
KB4561669	KB4486564	KB4019263
KB4557900	KB4486459	KB4018073
KB4556843	KB4483187	KB4015546
KB4550965	KB4480960	KB4012212
KB4541500	KB4474419	KB3205394
KB4537829	KB4471328	KB3197867
KB4537813	KB4467106	KB3192391
KB4536952	KB4462915	KB3172605
KB4534314	KB4457145	KB3161949

KB3159398	KB3061518	KB2937610
KB3156016	KB3060716	KB2931356
KB3155178	KB3059317	KB2923545
KB3150220	KB3046269	KB2912390
KB3139914	KB3046017	KB2911501
KB3139398	KB3045685	KB2900986
KB3138910	KB3042058	KB2894854
KB3138612	KB3037574	KB2894844
KB3133977	KB3035126	KB2893294
KB3127220	KB3031432	KB2888049
KB3126587	KB3030377	KB2884256
KB3124275	KB3023215	KB2882822
KB3122648	KB3022777	KB2871997
KB3115858	KB3021674	KB2868038
KB3110329	KB3020388	KB2864202
KB3109560	KB3019978	KB2862330
KB3109103	KB3011780	KB2862152
KB3108664	KB3010788	KB2861698
KB3108371	KB3004375	KB2857650
KB3101722	KB2992611	KB2852386
KB3097989	KB2991963	KB2847927
KB3093513	KB2984976	KB2840631
KB3092601	KB2984972	KB2834140
KB3086255	KB2978742	KB2830477
KB3080079	KB2977292	KB2813430
KB3075226	KB2973201	KB2813347
KB3075220	KB2973112	KB2807986
KB3074543	KB2972211	KB2786081
KB3072305	KB2972100	KB2770660
KB3071756	KB2968294	KB2758857
KB3067903	KB2943357	KB2742599

KB2736422

KB2654428

KB2564958

KB2731771

KB2653956

KB2560656

KB2729094

KB2639308

KB2552343

KB2727528

KB2631813

KB2533623

KB2705219

KB2621440

KB2533552

KB2698365

KB2620704

KB2532531

KB2690533

KB2604115

KB2506212

KB2685939

KB2585542

KB2491683

KB2670838

KB2579686

KB2479943

KB2667402

KB2574819

KB2467173

KB2656356

KB2565063